

TABLE 2

## EXAMPLE LOWER LIMITS OF QUANTITATION FOR SEMIVOLATILE ORGANICS

| Compound                     | Lower Limits of Quantitation <sup>a</sup> |   |
|------------------------------|---|---|
|                              | Ground water<br>(µg/L)                    | Low Soil/Sediment <sup>b</sup><br>(µg/kg) |
| Acenaphthene                 | 10  | 660                                       |
| Acenaphthylene               | 10  | 660                                       |
| Acetophenone                 | 10  | ND  |
| 2-Acetylaminofluorene        | 20  | ND  |
| 1-Acetyl-2-thiourea          | 1000                                      | ND  |
| 2-Aminoanthraquinone         | 20  | ND  |
| Aminoazobenzene              | 10  | ND  |
| 4-Aminobiphenyl              | 20  | ND  |
| Anilazine                    | 100                                       | ND  |
| o-Anisidine                  | 10  | ND  |
| Anthracene                   | 10  | 660                                       |
| Aramite                      | 20  | ND  |
| Azinphos-methyl              | 100                                       | ND  |
| Barban                       | 200                                       | ND  |
| Benz(a)anthracene            | 10  | 660                                       |
| Benzo(b)fluoranthene         | 10  | 660                                       |
| Benzo(k)fluoranthene         | 10  | 660                                       |
| Benzoic acid                 | 50  | 3300                                      |
| Benzo(g,h,i)perylene         | 10  | 660                                       |
| Benzo(a)pyrene               | 10  | 660                                       |
| p-Benzoquinone               | 10  | ND  |
| Benzyl alcohol               | 20  | 1300                                      |
| Bis(2-chloroethoxy)methane   | 10  | 660                                       |
| Bis(2-chloroethyl) ether     | 10  | 660                                       |
| Bis(2-chloroisopropyl) ether | 10  | 660                                       |
| 4-Bromophenyl phenyl ether   | 10  | 660                                       |
| Bromoxynil                   | 10  | ND  |
| Butyl benzyl phthalate       | 10  | 660                                       |
| Captafol                     | 20  | ND  |
| Captan                       | 50  | ND  |
| Carbaryl                     | 10  | ND  |
| Carbofuran                   | 10  | ND  |
| Carbophenothion              | 10  | ND  |
| Chlorfenvinphos              | 20  | ND  |
| 4-Chloroaniline              | 20  | 1300                                      |
| Chlorobenzilate              | 10  | ND  |
| 5-Chloro-2-methylaniline     | 10  | ND  |
| 4-Chloro-3-methylphenol      | 20  | 1300                                      |

TABLE 2  
(continued)

| Compound                               | Lower Limits of Quantitation <sup>a</sup> |   |
|--|---|---|
|  | Ground water<br>(µg/L)                    | Low Soil/Sediment <sup>b</sup><br>(µg/kg) |
| 3-(Chloromethyl)pyridine hydrochloride | 100                                       | ND  |
| 2-Chloronaphthalene                    | 10  | 660                                       |
| 2-Chlorophenol                         | 10  | 660                                       |
| 4-Chlorophenyl phenyl ether            | 10  | 660                                       |
| Chrysene                               | 10  | 660                                       |
| Coumaphos                              | 40  | ND  |
| p-Cresidine                            | 10  | ND  |
| Crotoxyphos                            | 20  | ND  |
| 2-Cyclohexyl-4,6-dinitrophenol         | 100                                       | ND  |
| Demeton-O                              | 10  | ND  |
| Demeton-S                              | 10  | ND  |
| Diallate (cis or trans)                | 10  | ND  |
| Diallate (trans or cis)                | 10  | ND  |
| 2,4-Diaminotoluene                     | 20  | ND  |
| Dibenz(a,j)acridine                    | 10  | ND  |
| Dibenz(a,h)anthracene                  | 10  | 660                                       |
| Dibenzofuran                           | 10  | 660                                       |
| Dibenzo(a,e)pyrene                     | 10  | ND  |
| Di-n-butyl phthalate                   | 10  | ND  |
| Dichlone                               | NA  | ND  |
| 1,2-Dichlorobenzene                    | 10  | 660                                       |
| 1,3-Dichlorobenzene                    | 10  | 660                                       |
| 1,4-Dichlorobenzene                    | 10  | 660                                       |
| 3,3'-Dichlorobenzidine                 | 20  | 1300                                      |
| 2,4-Dichlorophenol                     | 10  | 660                                       |
| 2,6-Dichlorophenol                     | 10  | ND  |
| Dichlorovos                            | 10  | ND  |
| Dicrotophos                            | 10  | ND  |
| Diethyl phthalate                      | 10  | 660                                       |
| Diethylstilbestrol                     | 20  | ND  |
| Diethyl sulfate                        | 100                                       | ND  |
| Dimethoate                             | 20  | ND  |
| 3,3'-Dimethoxybenzidine                | 100                                       | ND  |
| Dimethylaminoazobenzene                | 10  | ND  |
| 7,12-Dimethylbenz(a)anthracene         | 10  | ND  |
| 3,3'-Dimethylbenzidine                 | 10  | ND  |
| 2,4-Dimethylphenol                     | 10  | 660                                       |
| Dimethyl phthalate                     | 10  | 660                                       |
| 1,2-Dinitrobenzene                     | 40  | ND  |

TABLE 2  
(continued)

| Compound                    | Lower Limits of Quantitation <sup>a</sup> |   |
|-----------------------------|---|---|
|                             | Ground water<br>(µg/L)                    | Low Soil/Sediment <sup>b</sup><br>(µg/kg) |
| 1,3-Dinitrobenzene          | 20  | ND  |
| 1,4-Dinitrobenzene          | 40  | ND  |
| 4,6-Dinitro-2-methylphenol  | 50  | 3300                                      |
| 2,4-Dinitrophenol           | 50  | 3300                                      |
| 2,4-Dinitrotoluene          | 10  | 660                                       |
| 2,6-Dinitrotoluene          | 10  | 660                                       |
| Dinocap                     | 100                                       | ND  |
| Dinoseb                     | 20  | ND  |
| 5,5-Diphenylhydantoin       | 20  | ND  |
| Di-n-octyl phthalate        | 10  | 660                                       |
| Disulfoton                  | 10  | ND  |
| EPN                         | 10  | ND  |
| Ethion                      | 10  | ND  |
| Ethyl carbamate             | 50  | ND  |
| Bis(2-ethylhexyl) phthalate | 10  | 660                                       |
| Ethyl methanesulfonate      | 20  | ND  |
| Famphur                     | 20  | ND  |
| Fensulfothion               | 40  | ND  |
| Fenthion                    | 10  | ND  |
| Fluchloralin                | 20  | ND  |
| Fluoranthene                | 10  | 660                                       |
| Fluorene                    | 10  | 660                                       |
| Hexachlorobenzene           | 10  | 660                                       |
| Hexachlorobutadiene         | 10  | 660                                       |
| Hexachlorocyclopentadiene   | 10  | 660                                       |
| Hexachloroethane            | 10  | 660                                       |
| Hexachlorophene             | 50  | ND  |
| Hexachloropropene           | 10  | ND  |
| Hexamethylphosphoramide     | 20  | ND  |
| Indeno(1,2,3-cd)pyrene      | 10  | 660                                       |
| Isodrin                     | 20  | ND  |
| Isophorone                  | 10  | 660                                       |
| Isosafrole                  | 10  | ND  |
| Kepone                      | 20  | ND  |
| Leptophos                   | 10  | ND  |
| Malathion                   | 50  | ND  |
| Mestranol                   | 20  | ND  |
| Methapyrilene               | 100                                       | ND  |
| Methoxychlor                | 10  | ND  |

TABLE 2  
(continued)

| Compound                       | Lower Limits of Quantitation <sup>a</sup> |   |
|--------------------------------|---|---|
|                                | Ground water<br>(µg/L)                    | Low Soil/Sediment <sup>b</sup><br>(µg/kg) |
| 3-Methylcholanthrene           | 10  | ND  |
| Methyl methanesulfonate        | 10  | ND  |
| 2-Methylnaphthalene            | 10  | 660                                       |
| Methyl parathion               | 10  | ND  |
| 2-Methylphenol                 | 10  | 660                                       |
| 3-Methylphenol                 | 10  | ND  |
| 4-Methylphenol                 | 10  | 660                                       |
| Mevinphos                      | 10  | ND  |
| Mexacarbate                    | 20  | ND  |
| Mirex                          | 10  | ND  |
| Monocrotophos                  | 40  | ND  |
| Naled                          | 20  | ND  |
| Naphthalene                    | 10  | 660                                       |
| 1,4-Naphthoquinone             | 10  | ND  |
| 1-Naphthylamine                | 10  | ND  |
| 2-Naphthylamine                | 10  | ND  |
| Nicotine                       | 20  | ND  |
| 5-Nitroacenaphthene            | 10  | ND  |
| 2-Nitroaniline                 | 50  | 3300                                      |
| 3-Nitroaniline                 | 50  | 3300                                      |
| 4-Nitroaniline                 | 20  | ND  |
| 5-Nitro-o-anisidine            | 10  | ND  |
| Nitrobenzene                   | 10  | 660                                       |
| 4-Nitrobiphenyl                | 10  | ND  |
| Nitrofen                       | 20  | ND  |
| 2-Nitrophenol                  | 10  | 660                                       |
| 4-Nitrophenol                  | 50  | 3300                                      |
| 5-Nitro-o-toluidine            | 10  | ND  |
| 4-Nitroquinoline-1-oxide       | 40  | ND  |
| N-Nitrosodi-n-butylamine       | 10  | ND  |
| N-Nitrosodiethylamine          | 20  | ND  |
| N-Nitrosodiphenylamine         | 10  | 660                                       |
| N-Nitroso-di-n-propylamine     | 10  | 660                                       |
| N-Nitrosopiperidine            | 20  | ND  |
| N-Nitrosopyrrolidine           | 40  | ND  |
| Octamethyl pyrophosphoramidate | 200                                       | ND  |
| 4,4'-Oxydianiline              | 20  | ND  |
| Parathion                      | 10  | ND  |
| Pentachlorobenzene             | 10  | ND  |

TABLE 2  
(continued)

| Compound                          | Lower Limits of Quantitation <sup>a</sup> |   |
|-----------------------------------|---|---|
|                                   | Ground water<br>(µg/L)                    | Low Soil/Sediment <sup>b</sup><br>(µg/kg) |
| Pentachloronitrobenzene           | 20  | ND  |
| Pentachlorophenol                 | 50  | 3300                                      |
| Phenacetin                        | 20  | ND  |
| Phenanthrene                      | 10  | 660                                       |
| Phenobarbital                     | 10  | ND  |
| Phenol                            | 10  | 660                                       |
| 1,4-Phenylenediamine              | 10  | ND  |
| Phorate                           | 10  | ND  |
| Phosalone                         | 100                                       | ND  |
| Phosmet                           | 40  | ND  |
| Phosphamidon                      | 100                                       | ND  |
| Phthalic anhydride                | 100                                       | ND  |
| 2-Picoline                        | ND  | ND  |
| Piperonyl sulfoxide               | 100                                       | ND  |
| Pronamide                         | 10  | ND  |
| Propylthiouracil                  | 100                                       | ND  |
| Pyrene                            | 10  | 660                                       |
| Resorcinol                        | 100                                       | ND  |
| Safrole                           | 10  | ND  |
| Strychnine                        | 40  | ND  |
| Sulfallate                        | 10  | ND  |
| Terbufos                          | 20  | ND  |
| 1,2,4,5-Tetrachlorobenzene        | 10  | ND  |
| 2,3,4,6-Tetrachlorophenol         | 10  | ND  |
| Tetrachlorvinphos                 | 20  | ND  |
| Tetraethyl pyrophosphate          | 40  | ND  |
| Thionazine                        | 20  | ND  |
| Thiophenol (Benzenethiol)         | 20  | ND  |
| o-Toluidine                       | 10  | ND  |
| 1,2,4-Trichlorobenzene            | 10  | 660                                       |
| 2,4,5-Trichlorophenol             | 10  | 660                                       |
| 2,4,6-Trichlorophenol             | 10  | 660                                       |
| Trifluralin                       | 10  | ND  |
| 2,4,5-Trimethylaniline            | 10  | ND  |
| Trimethyl phosphate               | 10  | ND  |
| 1,3,5-Trinitrobenzene             | 10  | ND  |
| Tris(2,3-dibromopropyl) phosphate | 200                                       | ND  |
| Tri-p-tolyl phosphate(h)          | 10  | ND  |

TABLE 2  
(continued)

<sup>a</sup> Sample lower limits of quantitation are highly matrix-dependent and those listed here are provided for guidance and may not always be achievable.

<sup>b</sup> Lower limits of quantitation listed for soil/sediment are based on wet weight. When data are reported on a dry weight basis, the lower limits will be higher based on the % dry weight of each sample. These lower limits are based on a 30-g sample and gel permeation chromatography cleanup.

ND = Not Determined

NA = Not Applicable

Other Matrices

Factor<sup>c</sup>

High-concentration soil and sludges by ultrasonic extractor  
Non-water miscible waste

7.5  
75

<sup>c</sup>Lower limit of quantitation = (Lower limit of quantitation for low soil/sediment given above in Table 2) x (Factor)

TABLE 3

DFTPP KEY IONS AND ION ABUNDANCE CRITERIA<sup>a,b</sup>

| Mass | Ion Abundance Criteria          |
|------|---------------------------------|
| 51   | 10-80% of Base Peak             |
| 68   | < 2% of mass 69                 |
| 70   | < 2% of mass 69                 |
| 127  | 10-80% of Base Peak             |
| 197  | < 2% of mass 198                |
| 198  | Base peak, or > 50% of Mass 442 |
| 199  | 5-9% of mass 198                |
| 275  | 10-60% of Base Peak             |
| 365  | > 1% of mass 198                |
| 441  | present but < 24% of mass 442   |
| 442  | Base Peak, or > 50% of mass 198 |
| 443  | 15-24% of mass 442              |

<sup>a</sup> The majority of the data are taken from Reference 13 (Method 525.2).

<sup>b</sup> The criteria in this table are intended to be used as default criteria for quadrupole instrumentation if optimized manufacturer's operating conditions are not available. Alternate tuning criteria may be employed (e.g., CLP or Method 625), provided that method performance is not adversely affected. See Sec. 11.3.1

TABLE 4

RECOMMENDED MINIMUM RESPONSE FACTOR CRITERIA FOR INITIAL AND CONTINUING CALIBRATION VERIFICATION USING THE SUGGESTED IONS FROM TABLE 1

| Semivolatile Compounds        | Minimum Response Factor (RF) |
|-------------------------------|------------------------------|
| Benzaldehyde                  | 0.010                        |
| Phenol                        | 0.800                        |
| Bis(2-chloroethyl)ether       | 0.700                        |
| 2-Chlorophenol                | 0.800                        |
| 2-Methylphenol                | 0.700                        |
| 2,2'-Oxybis-(1-chloropropane) | 0.010                        |
| Acetophenone                  | 0.010                        |
| 4-Methylphenol                | 0.600                        |
| N-Nitroso-di-n-propylamine    | 0.500                        |
| Hexachloroethane              | 0.300                        |
| Nitrobenzene                  | 0.200                        |
| Isophorone                    | 0.400                        |
| 2-Nitrophenol                 | 0.100                        |
| 2,4-Dimethylphenol            | 0.200                        |
| Bis(2-chloroethoxy)methane    | 0.300                        |
| 2,4-Dichlorophenol            | 0.200                        |
| Naphthalene                   | 0.700                        |
| 4-Chloroaniline               | 0.010                        |
| Hexachlorobutadiene           | 0.010                        |
| Caprolactam                   | 0.010                        |
| 4-Chloro-3-methylphenol       | 0.200                        |
| 2-Methylnaphthalene           | 0.400                        |
| Hexachlorocyclopentadiene     | 0.050                        |
| 2,4,6-Trichlorophenol         | 0.200                        |
| 2,4,5-Trichlorophenol         | 0.200                        |
| 1,1'-Biphenyl                 | 0.010                        |
| 2-Chloronaphthalene           | 0.800                        |



TABLE 4  
(continued)

| Semivolatile Compounds      | Minimum Response Factor (RF) |
|-----------------------------|------------------------------|
| 2-Nitroaniline              | 0.010                        |
| Dimethyl phthalate          | 0.010                        |
| 2,6-Dinitrotoluene          | 0.200                        |
| Acenaphthylene              | 0.900                        |
| 3-Nitroaniline              | 0.010                        |
| Acenaphthene                | 0.900                        |
| 2,4-Dinitrophenol           | 0.010                        |
| 4-Nitrophenol               | 0.010                        |
| Dibenzofuran                | 0.800                        |
| 2,4-Dinitrotoluene          | 0.200                        |
| Diethyl phthalate           | 0.010                        |
| 1,2,4,5-Tetrachlorobenzene  | 0.010                        |
| 4-Chlorophenyl-phenyl ether | 0.400                        |
| Fluorene                    | 0.900                        |
| 4-Nitroaniline              | 0.010                        |
| 4,6-Dinitro-2-methylphenol  | 0.010                        |
| 4-Bromophenyl-phenyl ether  | 0.100                        |
| N-Nitrosodiphenylamine      | 0.010                        |
| Hexachlorobenzene           | 0.100                        |
| Atrazine                    | 0.010                        |
| Pentachlorophenol           | 0.050                        |
| Phenanthrene                | 0.700                        |
| Anthracene                  | 0.700                        |
| Carbazole                   | 0.010                        |
| Di-n-butyl phthalate        | 0.010                        |
| Fluoranthene                | 0.600                        |
| Pyrene                      | 0.600                        |
| Butyl benzyl phthalate      | 0.010                        |
| 3,3'-Dichlorobenzidine      | 0.010                        |
| Benzo(a)anthracene          | 0.800                        |

TABLE 4  
(continued)

| Semivolatile Compounds      | Minimum Response Factor (RF) |
|-----------------------------|------------------------------|
| Chrysene                    | 0.700                        |
| Bis-(2-ethylhexyl)phthalate | 0.010                        |
| Di-n-octyl phthalate        | 0.010                        |
| Benzo(b)fluoranthene        | 0.700                        |
| Benzo(k)fluoranthene        | 0.700                        |
| Benzo(a)pyrene              | 0.700                        |
| Indeno(1,2,3-cd)pyrene      | 0.500                        |
| Dibenz(a,h)anthracene       | 0.400                        |
| Benzo(g,h,i)perylene        | 0.500                        |
| 2,3,4,6-Tetrachlorophenol   | 0.010                        |

TABLE 5

SEMIVOLATILE INTERNAL STANDARDS WITH CORRESPONDING ANALYTES  
ASSIGNED FOR QUANTITATION

| 1,4-Dichlorobenzene-d <sub>4</sub> | Naphthalene-d <sub>8</sub>         | Acenaphthene-d <sub>10</sub> |
|------------------------------------|------------------------------------|------------------------------|
| Aniline                            | Acetophenone                       | Acenaphthene                 |
| Benzyl alcohol                     | Benzoic acid                       | Acenaphthylene               |
| Bis(2-chloroethyl) ether           | Bis(2-chloroethoxy)methane         | 1-Chloronaphthalene          |
| Bis(2-chloroisopropyl) ether       | 4-Chloroaniline                    | 2-Chloronaphthalene          |
| 2-Chlorophenol                     | 4-Chloro-3-methylphenol            | 4-Chlorophenyl phenyl ether  |
| 1,3-Dichlorobenzene                | 2,4-Dichlorophenol                 | Dibenzofuran                 |
| 1,4-Dichlorobenzene                | 2,6-Dichlorophenol                 | Diethyl phthalate            |
| 1,2-Dichlorobenzene                | α,α-Dimethyl-<br>phenethylamine    | Dimethyl phthalate           |
| Ethyl methanesulfonate             | 2,4-Dimethylphenol                 | 2,4-Dinitrophenol            |
| 2-Fluorophenol (surr)              | Hexachlorobutadiene                | 2,4-Dinitrotoluene           |
| Hexachloroethane                   | Isophorone                         | 2,6-Dinitrotoluene           |
| Methyl methanesulfonate            | 2-Methylnaphthalene                | Fluorene                     |
| 2-Methylphenol                     | Naphthalene                        | 2-Fluorobiphenyl (surr)      |
| 4-Methylphenol                     | Nitrobenzene                       | Hexachlorocyclopentadiene    |
| N-Nitrosodimethylamine             | Nitrobenzene-d <sub>8</sub> (surr) | 1-Naphthylamine              |
| N-Nitroso-di-n-propylamine         | 2-Nitrophenol                      | 2-Naphthylamine              |
| Phenol                             | N-Nitrosodi-n-butylamine           | 2-Nitroaniline               |
| Phenol-d <sub>6</sub> (surr)       | N-Nitrosopiperidine                | 3-Nitroaniline               |
| 2-Picoline                         | 1,2,4-Trichlorobenzene             | 4-Nitroaniline               |
|                                    |                                    | 4-Nitrophenol                |
|                                    |                                    | Pentachlorobenzene           |
|                                    |                                    | 1,2,4,5-Tetrachlorobenzene   |
|                                    |                                    | 2,3,4,6-Tetrachlorophenol    |
|                                    |                                    | 2,4,6-Tribromophenol (surr)  |
|                                    |                                    | 2,4,6-Trichlorophenol        |
|                                    |                                    | 2,4,5-Trichlorophenol        |

(surr) = surrogate

TABLE 5  
(continued)

| Phenanthrene-d <sub>10</sub> | Chrysene-d <sub>12</sub>         | Perylene-d <sub>12</sub>       |
|------------------------------|----------------------------------|--------------------------------|
| 4-Aminobiphenyl              | Benzidine                        | Benzo(b)fluoranthene           |
| Anthracene                   | Benzo(a)anthracene               | Benzo(k)fluoranthene           |
| 4-Bromophenyl phenyl ether   | Bis(2-ethylhexyl) phthalate      | Benzo(g,h,i)perylene           |
| Di-n-butyl phthalate         | Butyl benzyl phthalate           | Benzo(a)pyrene                 |
| 4,6-Dinitro-2-methylphenol   | Chrysene                         | Dibenz(a,j)acridine            |
| Diphenylamine                | 3,3'-Dichlorobenzidine           | Dibenz(a,h)anthracene          |
| Fluoranthene                 | p-Dimethyl aminoazobenzene       | 7,12-Dimethylbenz(a)anthracene |
| Hexachlorobenzene            | Pyrene                           | Di-n-octyl phthalate           |
| N-Nitrosodiphenylamine       | Terphenyl-d <sub>14</sub> (surr) | Indeno(1,2,3-cd) pyrene        |
| Pentachlorophenol            |                                  | 3-Methylcholanthrene           |
| Pentachloronitrobenzene      |                                  |                                |
| Phenacetin                   |                                  |                                |
| Phenanthrene                 |                                  |                                |
| Pronamide                    |                                  |                                |

(surr) = surrogate

TABLE 6

EXAMPLE SINGLE LABORATORY PERFORMANCE DATA<sup>a</sup>

| Compound                     | Test conc.<br>(µg/L) | × of 5<br>replicates<br>(µg/L) | % Recovery<br>of Avg. |
|------------------------------|----------------------|--------------------------------|-----------------------|
| Acenaphthene                 | 50                   | 46.7                           | 93.4                  |
| Acenaphthylene               | 50                   | 46.1                           | 92.2                  |
| Aniline                      | 50                   | 8.3                            | 16.7                  |
| Anthracene                   | 50                   | 48.4                           | 96.8                  |
| Benzoic acid                 | 50                   | 43.7                           | 87.4                  |
| Benz(a)anthracene            | 50                   | 49.6                           | 99.2                  |
| Benzo(b)fluoranthene         | 50                   | 49.8                           | 99.6                  |
| Benzo(k)fluoranthene         | 50                   | 50.6                           | 101                   |
| Benzo(a)pyrene               | 50                   | 47.7                           | 95.5                  |
| Benzo(g,h,i)perylene         | 50                   | 52.6                           | 105                   |
| Benzyl alcohol               | 50                   | 44.4                           | 88.8                  |
| Bis(2-chloroethyl) ether     | 50                   | 44.2                           | 88.4                  |
| Bis(2-chloroethoxy)methane   | 50                   | 46.6                           | 93.1                  |
| Bis(2-chloroisopropyl) ether | 50                   | 43.4                           | 86.8                  |
| Bis(2-ethylhexyl) phthalate  | 50                   | 50.2                           | 100                   |
| 4-Bromophenyl phenyl ether   | 50                   | 48.6                           | 97.2                  |
| Butyl benzyl phthalate       | 50                   | 49.6                           | 99.3                  |
| Carbazole                    | 50                   | 52.1                           | 104                   |
| 2-Chloroaniline              | 50                   | 38.9                           | 77.7                  |
| 4-Chloro-3-methylphenol      | 50                   | 47.3                           | 94.6                  |
| 2-Chloronaphthalene          | 50                   | 45.3                           | 90.8                  |
| 2-Chlorophenol               | 50                   | 43.1                           | 86.2                  |
| 4-Chlorophenyl phenyl ether  | 50                   | 47.3                           | 94.6                  |
| Chrysene                     | 50                   | 50.3                           | 101                   |
| Dibenzofuran                 | 50                   | 47.4                           | 94.7                  |
| Dibenz(a,h)anthracene        | 50                   | 51.6                           | 103                   |
| Di-n-butyl phthalate         | 50                   | 50.5                           | 101                   |
| 1,2-Dichlorobenzene          | 50                   | 35.8                           | 71.6                  |
| 1,3-Dichlorobenzene          | 50                   | 33.3                           | 66.7                  |
| 1,4-Dichlorobenzene          | 50                   | 34.4                           | 68.7                  |
| 3,3'-Dichlorobenzidine       | 50                   | 32.0                           | 64.0                  |
| 2,4-Dichlorophenol           | 50                   | 47.4                           | 94.8                  |
| Diethyl phthalate            | 50                   | 50.0                           | 99.9                  |
| Dimethyl phthalate           | 50                   | 48.5                           | 97.0                  |
| 2,4-Dimethylphenol           | 50                   | 31.2                           | 62.3                  |
| 4,6-Dinitro-2-methylphenol   | 50                   | 57.6                           | 115                   |
| 2,4-Dinitrophenol            | 50                   | 58.7                           | 117                   |
| 2,4-Dinitrotoluene           | 50                   | 51.3                           | 103                   |

TABLE 6  
(continued)

| Compound                  | Test conc.<br>(µg/L) | ✕•of 5<br>replicates<br>(µg/L) | % Recovery<br>of Avg. |
|---------------------------|----------------------|--------------------------------|-----------------------|
| 2,6-Dinitrotoluene        | 50                   | 50.2                           | 100                   |
| Di-n-octyl phthalate      | 50                   | 51.1                           | 102                   |
| Fluoranthene              | 50                   | 51.0                           | 102                   |
| Fluorene                  | 50                   | 48.5                           | 97.0                  |
| Hexachlorobenzene         | 50                   | 49.0                           | 97.9                  |
| Hexachlorobutadiene       | 50                   | 34.7                           | 69.5                  |
| Hexachlorocyclopentadiene | 50                   | 1.9                            | 3.8                   |
| Hexachloroethane          | 50                   | 29.9                           | 58.8                  |
| Indeno(1,2,3-cd)pyrene    | 50                   | 51.7                           | 103                   |
| Isophorone                | 50                   | 47.1                           | 94.3                  |
| 2-Methylnaphthalene       | 50                   | 44.7                           | 89.4                  |
| 2-Methylphenol            | 50                   | 41.7                           | 83.4                  |
| 4-Methylphenol            | 50                   | 42.6                           | 85.2                  |
| Naphthalene               | 50                   | 43.4                           | 86.8                  |
| 2-Nitroaniline            | 50                   | 48.4                           | 96.7                  |
| 3-Nitroaniline            | 50                   | 46.8                           | 93.6                  |
| 4-Nitroaniline            | 50                   | 56.1                           | 112                   |
| Nitrobenzene              | 50                   | 47.1                           | 94.1                  |
| 2-Nitrophenol             | 50                   | 47.3                           | 94.6                  |
| 4-Nitrophenol             | 50                   | 55.4                           | 111                   |
| N-Nitrosodiphenylamine    | 50                   | 46.7                           | 93.4                  |
| N-Nitroso-di-propylamine  | 50                   | 44.6                           | 89.3                  |
| Pentachlorophenol         | 50                   | 56.9                           | 114                   |
| Phenanthrene              | 50                   | 49.7                           | 99.4                  |
| Phenol                    | 50                   | 40.9                           | 81.8                  |
| Pyrene                    | 50                   | 49.2                           | 98.4                  |
| 1,2,4-Trichlorobenzene    | 50                   | 39.1                           | 78.2                  |
| 2,4,5-Trichlorophenol     | 50                   | 47.7                           | 95.4                  |
| 2,4,6-Trichlorophenol     | 50                   | 49.2                           | 98.4                  |

✕• = Average recovery for five initial demonstration of capability measurements, in µg/L

<sup>a</sup> Extraction using acidic pH only with a modified continuous liquid-liquid extractor with hydrophobic membrane according to Method 3520. These values are for guidance only. Appropriate derivation of acceptance criteria for similar extraction conditions may result in much different recovery ranges. See Method 8000 for information on developing and updating acceptance criteria for method performance.

TABLE 7

## EXTRACTION EFFICIENCY AND AQUEOUS STABILITY RESULTS

| Compound                               | Percent Recovery, Day 0 |     | Percent Recovery, Day 7 |     |
|--|-------------------------|-----|-------------------------|-----|
|  | Mean                    | RSD | Mean                    | RSD |
| 3-Amino-9-ethylcarbazole               | 80                      | 8   | 73                      | 3   |
| 4-Chloro-1,2-phenylenediamine          | 91                      | 1   | 108                     | 4   |
| 4-Chloro-1,3-phenylenediamine          | 84                      | 3   | 70                      | 3   |
| 1,2-Dibromo-3-chloropropane            | 97                      | 2   | 98                      | 5   |
| Dinoseb                                | 99                      | 3   | 97                      | 6   |
| Parathion                              | 100                     | 2   | 103                     | 4   |
| 4,4'-Methylenebis(N,N-dimethylaniline) | 108                     | 4   | 90                      | 4   |
| 5-Nitro-o-toluidine                    | 99                      | 10  | 93                      | 4   |
| 2-Picoline                             | 80                      | 4   | 83                      | 4   |
| Tetraethyl dithiopyrophosphate         | 92                      | 7   | 70                      | 1   |

Data taken from Reference 6.

TABLE 8

MEAN PERCENT RECOVERIES AND PERCENT RSD VALUES FOR SEMIVOLATILE ORGANIC FROM SPIKED CLAY SOIL AND TOPSOIL BY AUTOMATED SOXHLET EXTRACTION (METHOD 3541) WITH HEXANE-ACETONE (1:1)<sup>a</sup>

| Compound                                   | Clay Soil     |     | Topsoil       |     |
|--|---------------|-----|---------------|-----|
|  | Mean Recovery | RSD | Mean Recovery | RSD |
| 1,3-Dichlorobenzene                        | 0             | --  | 0             | --  |
| 1,2-Dichlorobenzene                        | 0             | --  | 0             | --  |
| Nitrobenzene                               | 0             | --  | 0             | --  |
| Benzal chloride                            | 0             | --  | 0             | --  |
| Benzotrichloride                           | 0             | --  | 0             | --  |
| 4-Chloro-2-nitrotoluene                    | 0             | --  | 0             | --  |
| Hexachlorocyclopentadiene                  | 4.1           | 15  | 7.8           | 23  |
| 2,4-Dichloronitrobenzene                   | 35.2          | 7.6 | 21.2          | 15  |
| 3,4-Dichloronitrobenzene                   | 34.9          | 15  | 20.4          | 11  |
| Pentachlorobenzene                         | 13.7          | 7.3 | 14.8          | 13  |
| 2,3,4,5-Tetrachloronitrobenzene            | 55.9          | 6.7 | 50.4          | 6.0 |
| Benefin                                    | 62.6          | 4.8 | 62.7          | 2.9 |
| alpha-BHC                                  | 58.2          | 7.3 | 54.8          | 4.8 |
| Hexachlorobenzene                          | 26.9          | 13  | 25.1          | 5.7 |
| delta-BHC                                  | 95.8          | 4.6 | 99.2          | 1.3 |
| Heptachlor                                 | 46.9          | 9.2 | 49.1          | 6.3 |
| Aldrin                                     | 97.7          | 12  | 102           | 7.4 |
| Isopropalin                                | 102           | 4.3 | 105           | 2.3 |
| Heptachlor epoxide                         | 90.4          | 4.4 | 93.6          | 2.4 |
| trans-Chlordane                            | 90.1          | 4.5 | 95.0          | 2.3 |
| Endosulfan I                               | 96.3          | 4.4 | 101           | 2.2 |
| Dieldrin                                   | 129           | 4.7 | 104           | 1.9 |
| 2,5-Dichlorophenyl-4-nitrophenyl ether     | 110           | 4.1 | 112           | 2.1 |
| Endrin                                     | 102           | 4.5 | 106           | 3.7 |
| Endosulfan II                              | 104           | 4.1 | 105           | 0.4 |
| p,p'-DDT                                   | 134           | 2.1 | 111           | 2.0 |
| 2,3,6-Trichlorophenyl-4'-nitrophenyl ether | 110           | 4.8 | 110           | 2.8 |
| 2,3,4-Trichlorophenyl-4'-nitrophenyl ether | 112           | 4.4 | 112           | 3.3 |
| Mirex                                      | 104           | 5.3 | 108           | 2.2 |

<sup>a</sup> The operating conditions for the Soxtec apparatus were as follows: immersion time 45 min; extraction time 45 min; the sample size was 10 g; the spiking concentration was 500 ng/g, except for the surrogate compounds at 1000 ng/g, 2,5-Dichlorophenyl-4-nitrophenyl ether, 2,3,6-Trichlorophenyl-4-nitrophenyl ether, and 2,3,4-Trichlorophenyl-4-nitrophenyl ether at 1500 ng/g, Nitrobenzene at 2000 ng/g, and 1,3-Dichlorobenzene and 1,2-Dichlorobenzene at 5000 ng/g.



TABLE 9

SINGLE LABORATORY ACCURACY AND PRECISION DATA FOR THE EXTRACTION  
OF SEMIVOLATILE ORGANICS FROM SPIKED CLAY BY  
AUTOMATED SOXHLET (METHOD 3541)<sup>a</sup>

| Compound                    | Mean Recovery | RSD |
|-----------------------------|---------------|-----|
| Phenol                      | 47.8          | 5.6 |
| Bis(2-chloroethyl)ether     | 25.4          | 13  |
| 2-Chlorophenol              | 42.7          | 4.3 |
| Benzyl alcohol              | 55.9          | 7.2 |
| 2-Methylphenol              | 17.6          | 6.6 |
| Bis(2-chloroisopropyl)ether | 15.0          | 15  |
| 4-Methylphenol              | 23.4          | 6.7 |
| N-Nitroso-di-n-propylamine  | 41.4          | 6.2 |
| Nitrobenzene                | 28.2          | 7.7 |
| Isophorone                  | 56.1          | 4.2 |
| 2-Nitrophenol               | 36.0          | 6.5 |
| 2,4-Dimethylphenol          | 50.1          | 5.7 |
| Benzoic acid                | 40.6          | 7.7 |
| Bis(2-chloroethoxy)methane  | 44.1          | 3.0 |
| 2,4-Dichlorophenol          | 55.6          | 4.6 |
| 1,2,4-Trichlorobenzene      | 18.1          | 31  |
| Naphthalene                 | 26.2          | 15  |
| 4-Chloroaniline             | 55.7          | 12  |
| 4-Chloro-3-methylphenol     | 65.1          | 5.1 |
| 2-Methylnaphthalene         | 47.0          | 8.6 |
| Hexachlorocyclopentadiene   | 19.3          | 19  |
| 2,4,6-Trichlorophenol       | 70.2          | 6.3 |
| 2,4,5-Trichlorophenol       | 26.8          | 2.9 |
| 2-Chloronaphthalene         | 61.2          | 6.0 |
| 2-Nitroaniline              | 73.8          | 6.0 |
| Dimethyl phthalate          | 74.6          | 5.2 |
| Acenaphthylene              | 71.6          | 5.7 |
| 3-Nitroaniline              | 77.6          | 5.3 |
| Acenaphthene                | 79.2          | 4.0 |
| 2,4-Dinitrophenol           | 91.9          | 8.9 |
| 4-Nitrophenol               | 62.9          | 16  |
| Dibenzofuran                | 82.1          | 5.9 |
| 2,4-Dinitrotoluene          | 84.2          | 5.4 |
| 2,6-Dinitrotoluene          | 68.3          | 5.8 |

| Compound                    | Mean Recovery | RSD |
|-----------------------------|---------------|-----|
| Diethyl phthalate           | 74.9          | 5.4 |
| 4-Chlorophenyl-phenyl ether | 67.2          | 3.2 |
| Fluorene                    | 82.1          | 3.4 |
| 4-Nitroaniline              | 79.0          | 7.9 |
| 4,6-Dinitro-2-methylphenol  | 63.4          | 6.8 |
| N-Nitrosodiphenylamine      | 77.0          | 3.4 |
| 4-Bromophenyl-phenyl ether  | 62.4          | 3.0 |
| Hexachlorobenzene           | 72.6          | 3.7 |
| Pentachlorophenol           | 62.7          | 6.1 |
| Phenanthrene                | 83.9          | 5.4 |
| Anthracene                  | 96.3          | 3.9 |
| Di-n-butyl phthalate        | 78.3          | 40  |
| Fluoranthene                | 87.7          | 6.9 |
| Pyrene                      | 102           | 0.8 |
| Butyl benzyl phthalate      | 66.3          | 5.2 |
| 3,3'-Dichlorobenzidine      | 25.2          | 11  |
| Benzo(a)anthracene          | 73.4          | 3.8 |
| Bis(2-ethylhexyl) phthalate | 77.2          | 4.8 |
| Chrysene                    | 76.2          | 4.4 |
| Di-n-octyl phthalate        | 83.1          | 4.8 |
| Benzo(b)fluoranthene        | 82.7          | 5.0 |
| Benzo(k)fluoranthene        | 71.7          | 4.1 |
| Benzo(a)pyrene              | 71.7          | 4.1 |
| Indeno(1,2,3-cd)pyrene      | 72.2          | 4.3 |
| Dibenz(a,h)anthracene       | 66.7          | 6.3 |
| Benzo(g,h,i)perylene        | 63.9          | 8.0 |
| 1,2-Dichlorobenzene         | 0             | --  |
| 1,3-Dichlorobenzene         | 0             | --  |
| 1,4-Dichlorobenzene         | 0             | --  |
| Hexachloroethane            | 0             | --  |
| Hexachlorobutadiene         | 0             | --  |

<sup>a</sup> Number of determinations was three. The operating conditions for the Soxtec apparatus were as follows: immersion time 45 min; extraction time 45 min; the sample size was 10 g clay soil; the spike concentration was 6 mg/kg per compound. The sample was allowed to equilibrate 1 hour after spiking.

Data taken from Reference 7.

TABLE 10

PRECISION AND BIAS VALUES FOR METHOD 3542<sup>1</sup>

| Compound                    | Mean Recovery | Standard Deviation | % RSD |
|-----------------------------|---------------|--------------------|-------|
| 2-Fluorophenol              | 74.6          | 28.6               | 38.3  |
| Phenol-d <sub>5</sub>       | 77.8          | 27.7               | 35.6  |
| Nitrobenzene-d <sub>5</sub> | 65.6          | 32.5               | 49.6  |
| 2-Fluorobiphenyl            | 75.9          | 30.3               | 39.9  |
| 2,4,6-Tribromophenol        | 67.0          | 34.0               | 50.7  |
| Terphenyl-d <sub>14</sub>   | 78.6          | 32.4               | 41.3  |

<sup>1</sup> The surrogate values shown in Table 10 represent mean recoveries for surrogates in all Method 0010 matrices in a field dynamic spiking study.

TABLE 11

PRESSURIZED FLUID EXTRACTION (METHOD 3545) RECOVERY VALUES  
AS PERCENT OF SOXTEC™

| Compound                    | Clay   |        |        | Loam   |        |        | Sand   |        |        | Mean<br>Rec. |
|-----------------------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------------|
|                             | Low    | Mid    | High   | Low    | Mid    | High   | Low    | Mid    | High   |              |
| Phenol                      | 93.3   | 78.7   | 135.9  | 73.9   | 82.8   | 124.6  | 108.8  | 130.6  | 89.7   | 102.0        |
| Bis(2-chloroethyl) ether    | 102.1  | 85.1   | 109.1  | 96.0   | 88.0   | 103.6  | 122.3  | 119.9  | 90.8   | 101.9        |
| 2-Chlorophenol              | 100.8  | 82.6   | 115.0  | 93.8   | 88.9   | 111.1  | 115.0  | 115.3  | 91.9   | 101.6        |
| 1,3-Dichlorobenzene         | 127.7  | 129.7  | 110.0  | *364.2 | 129.9  | 119.0  | *241.3 | *163.7 | 107.1  | 120.6        |
| 1,4-Dichlorobenzene         | 127.9  | 127.0  | 110.5  | *365.9 | 127.8  | 116.4  | *309.6 | *164.1 | 105.8  | 119.2        |
| 1,2-Dichlorobenzene         | 116.8  | 115.8  | 101.3  | *159.2 | 113.4  | 105.5  | *189.3 | 134.0  | 100.4  | 112.5        |
| 2-Methylphenol              | 98.9   | 82.1   | 119.7  | 87.6   | 89.4   | 111.0  | 133.2  | 128.0  | 92.1   | 104.7        |
| Bis(2-chloroisopropyl)ether | 109.4  | 71.5   | 108.0  | 81.8   | 81.0   | 88.6   | 118.1  | 148.3  | 94.8   | 100.2        |
| o-Toluidine                 | 100.0  | 89.7   | 117.2  | 100.0  | *152.5 | 120.3  | 100.0  | *199.5 | 102.7  | 110.3        |
| N-Nitroso-di-n-propylamine  | 103.0  | 79.1   | 107.7  | 83.9   | 88.1   | 96.2   | 109.9  | 123.3  | 91.4   | 98.1         |
| Hexachloroethane            | 97.1   | 125.1  | 111.0  | *245.4 | 117.1  | 128.1  | *566.7 | 147.9  | 103.7  | 118.6        |
| Nitrobenzene                | 104.8  | 82.4   | 106.6  | 86.8   | 84.6   | 101.7  | 119.7  | 122.1  | 93.3   | 100.2        |
| Isophorone                  | 100.0  | 86.4   | 98.2   | 87.1   | 87.5   | 109.7  | 135.5  | 118.4  | 92.7   | 101.7        |
| 2,4-Dimethylphenol          | 100.0  | 104.5  | 140.0  | 100.0  | 114.4  | 123.1  | 100.0  | *180.6 | 96.3   | 109.8        |
| 2-Nitrophenol               | 80.7   | 80.5   | 107.9  | 91.4   | 86.7   | 103.2  | 122.1  | 107.1  | 87.0   | 96.3         |
| Bis(chloroethoxy)methane    | 94.4   | 80.6   | 94.7   | 86.5   | 84.4   | 99.6   | 130.6  | 110.7  | 93.2   | 97.2         |
| 2,4-Dichlorophenol          | 88.9   | 87.8   | 111.4  | 85.9   | 87.6   | 103.5  | 123.3  | 107.0  | 92.1   | 98.6         |
| 1,2,4-Trichlorobenzene      | 98.0   | 97.8   | 98.8   | 123.0  | 93.7   | 94.5   | 137.0  | 99.4   | 95.3   | 104.2        |
| Naphthalene                 | 101.7  | 97.2   | 123.6  | 113.2  | 102.9  | 129.5  | *174.5 | 114.0  | 89.8   | 106.1        |
| 4-Chloroaniline             | 100.0  | *150.2 | *162.4 | 100.0  | 125.5  | *263.6 | 100.0  | *250.8 | 114.9  | 108.1        |
| Hexachlorobutadiene         | 101.1  | 98.7   | 102.2  | 124.1  | 90.3   | 98.0   | 134.9  | 96.1   | 96.8   | 104.7        |
| 4-Chloro-3-methylphenol     | 90.4   | 80.2   | 114.7  | 79.0   | 85.2   | 109.8  | 131.6  | 116.2  | 90.1   | 99.7         |
| 2-Methylnaphthalene         | 93.2   | 89.9   | 94.6   | 104.1  | 92.2   | 105.9  | 146.2  | 99.1   | 93.3   | 102.1        |
| Hexachlorocyclopentadiene   | 100.0  | 100.0  | 0.0    | 100.0  | 100.0  | 6.8    | 100.0  | 100.0  | *238.3 | 75.8         |
| 2,4,6-Trichlorophenol       | 94.6   | 90.0   | 112.0  | 84.2   | 91.2   | 103.6  | 101.6  | 95.9   | 89.8   | 95.9         |
| 2,4,5-Trichlorophenol       | 84.4   | 91.9   | 109.6  | 96.1   | 80.7   | 103.6  | 108.9  | 83.9   | 87.9   | 94.1         |
| 2-Chloronaphthalene         | 100.0  | 91.3   | 93.6   | 97.6   | 93.4   | 98.3   | 106.8  | 93.0   | 92.0   | 96.2         |
| 2-Nitroaniline              | 90.0   | 83.4   | 97.4   | 71.3   | 88.4   | 89.9   | 112.1  | 113.3  | 87.7   | 92.6         |
| 2,6-Dinitrotoluene          | 83.1   | 90.6   | 91.6   | 86.4   | 90.6   | 90.3   | 104.3  | 84.7   | 90.9   | 90.3         |
| Acenaphthylene              | 104.9  | 95.9   | 100.5  | 99.0   | 97.9   | 108.8  | 118.5  | 97.8   | 92.0   | 101.7        |
| 3-Nitroaniline              | *224.0 | 115.6  | 97.6   | 100.0  | 111.8  | 107.8  | 0.0    | 111.7  | 99.0   | 92.9         |
| Acenaphthene                | 102.1  | 92.6   | 97.6   | 97.2   | 96.9   | 104.4  | 114.2  | 92.0   | 89.0   | 98.4         |
| 4-Nitrophenol               | 0.0    | 93.2   | 121.5  | 18.1   | 87.1   | 116.6  | 69.1   | 90.5   | 84.5   | 75.6         |
| 2,4-Dinitrotoluene          | 73.9   | 91.9   | 100.2  | 84.7   | 93.8   | 98.9   | 100.9  | 84.3   | 87.3   | 90.7         |

TABLE 11  
(continued)

| Compound                    | Clay  |        |       | Loam  |        |        | Sand   |        |       | Mean Rec. |
|-----------------------------|-------|--------|-------|-------|--------|--------|--------|--------|-------|-----------|
|                             | Low   | Mid    | High  | Low   | Mid    | High   | Low    | Mid    | High  |           |
| Dibenzofuran                | 89.5  | 91.7   | 109.3 | 98.5  | 92.2   | 111.4  | 113.8  | 92.7   | 90.4  | 98.8      |
| 4-Chlorophenyl phenyl ether | 83.0  | 94.5   | 98.7  | 95.7  | 94.3   | 94.2   | 111.4  | 87.7   | 90.3  | 94.4      |
| Fluorene                    | 85.2  | 94.9   | 89.2  | 102.0 | 95.5   | 93.8   | 121.3  | 85.7   | 90.9  | 95.4      |
| 4-Nitroaniline              | 77.8  | 114.8  | 94.5  | 129.6 | 103.6  | 95.4   | *154.1 | 89.3   | 87.5  | 99.1      |
| N-Nitrosodiphenylamine      | 82.6  | 96.7   | 93.8  | 92.9  | 93.4   | 116.4  | 97.5   | 110.9  | 86.7  | 96.8      |
| 4-Bromophenyl phenyl ether  | 85.6  | 92.9   | 92.8  | 91.1  | 107.6  | 89.4   | 118.0  | 97.5   | 87.1  | 95.8      |
| Hexachlorobenzene           | 95.4  | 91.7   | 92.3  | 95.4  | 93.6   | 83.7   | 106.8  | 94.3   | 90.0  | 93.7      |
| Pentachlorophenol           | 68.2  | 85.9   | 107.7 | 53.2  | 89.8   | 88.1   | 96.6   | 59.8   | 81.3  | 81.2      |
| Phenanthrene                | 92.1  | 93.7   | 93.3  | 100.0 | 97.8   | 113.3  | 124.4  | 101.0  | 89.9  | 100.6     |
| Anthracene                  | 101.6 | 95.0   | 93.5  | 92.5  | 101.8  | 118.4  | 123.0  | 94.5   | 90.6  | 101.2     |
| Carbazole                   | 94.4  | 99.3   | 96.6  | 105.5 | 96.7   | 111.4  | 115.7  | 83.2   | 88.9  | 99.1      |
| Fluoranthene                | 109.9 | 101.4  | 94.3  | 111.6 | 96.6   | 109.6  | 123.2  | 85.4   | 92.7  | 102.7     |
| Pyrene                      | 106.5 | 105.8  | 107.6 | 116.7 | 90.7   | 127.5  | 103.4  | 95.5   | 93.2  | 105.2     |
| 3,3'-Dichlorobenzidine      | 100.0 | *492.3 | 131.4 | 100.0 | *217.6 | *167.6 | 100.0  | *748.8 | 100.0 | 116.5     |
| Benzo(a)anthracene          | 98.1  | 107.0  | 98.4  | 119.3 | 98.6   | 104.0  | 105.0  | 93.4   | 89.3  | 101.5     |
| Chrysene                    | 100.0 | 108.5  | 100.2 | 116.8 | 93.0   | 117.0  | 106.7  | 93.6   | 90.2  | 102.9     |
| Benzo(b)fluoranthene        | 106.6 | 109.9  | 75.6  | 121.7 | 100.7  | 93.9   | 106.9  | 81.9   | 93.6  | 99.0      |
| Benzo(k)fluoranthene        | 102.4 | 105.2  | 88.4  | 125.5 | 99.4   | 95.1   | 144.7  | 89.2   | 78.1  | 103.1     |
| Benzo(a)pyrene              | 107.9 | 105.5  | 80.8  | 122.3 | 97.7   | 104.6  | 101.7  | 86.2   | 92.0  | 99.9      |
| Indeno(1,2,3-cd)pyrene      | 95.1  | 105.7  | 93.8  | 126.0 | 105.2  | 90.4   | 133.6  | 82.6   | 91.9  | 102.7     |
| Dibenz(a,h)anthracene       | 85.0  | 102.6  | 82.0  | 118.8 | 100.7  | 91.9   | 142.3  | 71.0   | 93.1  | 98.6      |
| Benzo(g,h,i)perylene        | 98.0  | 0.0    | 81.2  | 0.0   | 33.6   | 78.6   | 128.7  | 83.0   | 94.2  | 66.4      |
| Mean                        | 95.1  | 94.3   | 101.0 | 95.5  | 96.5   | 104.1  | 113.0  | 100.9  | 92.5  |           |

\* Values greater than 150% were not used to determine the averages, but the 0% values were used.

TABLE 12

SINGLE LABORATORY ACCURACY AND PRECISION FOR THE EXTRACTION OF PAHs  
FROM A CERTIFIED REFERENCE SEDIMENT EC-1, USING METHOD 3561  
(SFE - SOLID TRAP)

| Compound               | Certified Value<br>(mg/kg) | SFE Value <sup>a</sup><br>(mg/kg) | Percent of<br>Certified Value | SFE<br>RSD |
|------------------------|----------------------------|-----------------------------------|-------------------------------|------------|
| Naphthalene            | (27.9) <sup>b</sup>        | 41.3 ± 3.6                        | (148)                         | 8.7        |
| Acenaphthylene         | (0.8)                      | 0.9 ± 0.1                         | (112)                         | 11.1       |
| Acenaphthene           | (0.2)                      | 0.2 ± 0.01                        | (100)                         | 0.05       |
| Fluorene               | (15.3)                     | 15.6 ± 1.8                        | (102)                         | 11.5       |
| Phenanthrene           | 15.8 ± 1.2                 | 16.1 ± 1.8                        | 102                           | 11.2       |
| Anthracene             | (1.3)                      | 1.1 ± 0.2                         | (88)                          | 18.2       |
| Fluoranthene           | 23.2 ± 2.0                 | 24.1 ± 2.1                        | 104                           | 8.7        |
| Pyrene                 | 16.7 ± 2.0                 | 17.2 ± 1.9                        | 103                           | 11.0       |
| Benz(a)anthracene      | 8.7 ± 0.8                  | 8.8 ± 1.0                         | 101                           | 11.4       |
| Chrysene               | (9.2)                      | 7.9 ± 0.9                         | (86)                          | 11.4       |
| Benzo(b)fluoranthene   | 7.9 ± 0.9                  | 8.5 ± 1.1                         | 108                           | 12.9       |
| Benzo(k)fluoranthene   | 4.4 ± 0.5                  | 4.1 ± 0.5                         | 91                            | 12.2       |
| Benzo(a)pyrene         | 5.3 ± 0.7                  | 5.1 ± 0.6                         | 96                            | 11.8       |
| Indeno(1,2,3-cd)pyrene | 5.7 ± 0.6                  | 5.2 ± 0.6                         | 91                            | 11.5       |
| Benzo(g,h,i)perylene   | 4.9 ± 0.7                  | 4.3 ± 0.5                         | 88                            | 11.6       |
| Dibenz(a,h)anthracene  | (1.3)                      | 1.1 ± 0.2                         | (85)                          | 18.2       |

<sup>a</sup> Relative standard deviations for the SFE values are based on six replicate extractions.

<sup>b</sup> Values in parentheses were obtained from, or compared to, Soxhlet extraction results which were not certified.

Data are taken from Reference 10.

TABLE 13

SINGLE LABORATORY ACCURACY AND PRECISION FOR THE EXTRACTION OF PAHs  
FROM A CERTIFIED REFERENCE SEDIMENT HS-3, USING METHOD 3561  
(SFE - SOLID TRAP)

| Compound               | Certified Value<br>(mg/kg) | SFE Value <sup>a</sup><br>(mg/kg) | Percent of<br>Certified Value | SFE<br>RSD |
|------------------------|----------------------------|-----------------------------------|-------------------------------|------------|
| Naphthalene            | 9.0 ± 0.7                  | 7.4 ± 0.6                         | 82                            | 8.1        |
| Acenaphthylene         | 0.3 ± 0.1                  | 0.4 ± 0.1                         | 133                           | 25.0       |
| Acenaphthene           | 4.5 ± 1.5                  | 3.3 ± 0.3                         | 73                            | 9.0        |
| Fluorene               | 13.6 ± 3.1                 | 10.4 ± 1.3                        | 77                            | 12.5       |
| Phenanthrene           | 85.0 ± 20.0                | 86.2 ± 9.5                        | 101                           | 11.0       |
| Anthracene             | 13.4 ± 0.5                 | 12.1 ± 1.5                        | 90                            | 12.4       |
| Fluoranthene           | 60.0 ± 9.0                 | 54.0 ± 6.1                        | 90                            | 11.3       |
| Pyrene                 | 39.0 ± 9.0                 | 32.7 ± 3.7                        | 84                            | 11.3       |
| Benz(a)anthracene      | 14.6 ± 2.0                 | 12.1 ± 1.3                        | 83                            | 10.7       |
| Chrysene               | 14.1 ± 2.0                 | 12.0 ± 1.3                        | 85                            | 10.8       |
| Benzo(b)fluoranthene   | 7.7 ± 1.2                  | 8.4 ± 0.9                         | 109                           | 10.7       |
| Benzo(k)fluoranthene   | 2.8 ± 2.0                  | 3.2 ± 0.5                         | 114                           | 15.6       |
| Benzo(a)pyrene         | 7.4 ± 3.6                  | 6.6 ± 0.8                         | 89                            | 12.1       |
| Indeno(1,2,3-cd)pyrene | 5.0 ± 2.0                  | 4.5 ± 0.6                         | 90                            | 13.3       |
| Benzo(g,h,i)perylene   | 5.4 ± 1.3                  | 4.4 ± 0.6                         | 82                            | 13.6       |
| Dibenz(a,h)anthracene  | 1.3 ± 0.5                  | 1.1 ± 0.3                         | 85                            | 27.3       |

<sup>a</sup> Relative standard deviations for the SFE values are based on three replicate extractions.

Data are taken from Reference 10.

TABLE 14

SINGLE LABORATORY ACCURACY AND PRECISION FOR THE EXTRACTION OF PAHs  
FROM A CERTIFIED REFERENCE SOIL SRS103-100, USING METHOD 3561  
(SFE - LIQUID TRAP)

| Compound                                       | Certified Value<br>(mg/kg) | SFE Value <sup>a</sup><br>(mg/kg) | Percent of<br>Certified Value | SFE<br>RSD |
|--|----------------------------|-----------------------------------|-------------------------------|------------|
| Naphthalene                                    | 32.4 ± 8.2                 | 29.55                             | 91                            | 10.5       |
| 2-Methylnaphthalene                            | 62.1 ± 11.5                | 76.13                             | 122                           | 2.0        |
| Acenaphthene                                   | 632 ± 105                  | 577.28                            | 91                            | 2.9        |
| Dibenzofuran                                   | 307 ± 49                   | 302.25                            | 98                            | 4.1        |
| Fluorene                                       | 492 ± 78                   | 427.15                            | 87                            | 3.0        |
| Phenanthrene                                   | 1618 ± 340                 | 1278.03                           | 79                            | 3.4        |
| Anthracene                                     | 422 ± 49                   | 400.80                            | 95                            | 2.6        |
| Fluoranthene                                   | 1280 ± 220                 | 1019.13                           | 80                            | 4.5        |
| Pyrene   | 1033 ± 285                 | 911.82                            | 88                            | 3.1        |
| Benz(a)anthracene                              | 252 ± 8                    | 225.50                            | 89                            | 4.8        |
| Chrysene                                       | 297 ± 26                   | 283.00                            | 95                            | 3.8        |
| Benzo(a)pyrene                                 | 97.2 ± 17.1                | 58.28                             | 60                            | 6.5        |
| Benzo(b)fluoranthene +<br>Benzo(k)fluoranthene | 153 ± 22                   | 130.88                            | 86                            | 10.7       |

<sup>a</sup> Relative standard deviations for the SFE values are based on four replicate extractions.

Data are taken from Reference 11.



TABLE 15

SINGLE LABORATORY RECOVERY DATA FOR SOLID-PHASE EXTRACTION (METHOD  
3535) OF BASE/NEUTRAL/ACID EXTRACTABLES FROM SPIKED TCLP BUFFERS  
LOW SPIKE LEVEL

| Analyte               | Spike<br>Level<br>(µg/L) | Buffer 1 (pH = 2.886) |     | Buffer 2 (pH = 4.937) |     |
|-----------------------|--------------------------|-----------------------|-----|-----------------------|-----|
|                       |                          | Recovery (%)          | RSD | Recovery (%)          | RSD |
| 1,4-Dichlorobenzene   | 3,750                    | 63                    | 10  | 63                    | 9   |
| Hexachloroethane      | 1,500                    | 55                    | 6   | 77                    | 4   |
| Nitrobenzene          | 1,000                    | 82                    | 10  | 100                   | 5   |
| Hexachlorobutadiene   | 250                      | 65                    | 3   | 56                    | 4   |
| 2,4-Dinitrotoluene    | 65                       | 89                    | 4   | 101                   | 5   |
| Hexachlorobenzene     | 65                       | 98                    | 5   | 95                    | 6   |
| o-Cresol              | 100,000                  | 83                    | 10  | 85                    | 5   |
| m-Cresol*             | 100,000                  | 86                    | 8   | 85                    | 3   |
| p-Cresol*             | 100,000                  | *                     | *   | *                     | *   |
| 2,4,6-Trichlorophenol | 1,000                    | 84                    | 12  | 95                    | 12  |
| 2,4,5-Trichlorophenol | 200,000                  | 83                    | 11  | 88                    | 3   |
| Pentachlorophenol     | 50,000                   | 82                    | 9   | 78                    | 9   |

Results from seven replicate spiked buffer samples.

\* In this study, m-cresol and p-cresol co-eluted and were quantitated as a mixture of both isomers.

Data from Reference 12.

TABLE 16

SINGLE LABORATORY RECOVERY DATA FOR SOLID-PHASE EXTRACTION (METHOD  
3535) OF BASE/NEUTRAL/ACID EXTRACTABLES FROM SPIKED TCLP BUFFERS  
HIGH SPIKE LEVEL

| Analyte               | Spike<br>Level<br>(µg/L) | Buffer 1 (pH = 2.886) |     | Buffer 2 (pH = 4.937) |     |
|-----------------------|--------------------------|-----------------------|-----|-----------------------|-----|
|                       |                          | Recovery (%)          | RSD | Recovery (%)          | RSD |
| 1,4-Dichlorobenzene   | 15,000                   | 63                    | 10  | 63                    | 9   |
| Hexachloroethane      | 6,000                    | 54                    | 7   | 46                    | 7   |
| Nitrobenzene          | 4,000                    | 81                    | 4   | 81                    | 13  |
| Hexachlorobutadiene   | 1,000                    | 81                    | 5   | 70                    | 11  |
| 2,4-Dinitrotoluene    | 260                      | 99                    | 8   | 98                    | 3   |
| Hexachlorobenzene     | 260                      | 89                    | 8   | 91                    | 9   |
| o-Cresol*             | 400,000                  | 92                    | 15  | 90                    | 4   |
| m-Cresol*             | 400,000                  | 95                    | 8   | 82                    | 6   |
| p-Cresol*             | 400,000                  | 82                    | 14  | 84                    | 7   |
| 2,4,6-Trichlorophenol | 4,000                    | 93                    | 12  | 104                   | 12  |
| 2,4,5-Trichlorophenol | 800,000                  | 93                    | 14  | 97                    | 23  |
| Pentachlorophenol     | 200,000                  | 84                    | 9   | 73                    | 8   |

Results from seven replicate spiked buffer samples.

\* In this study, recoveries of these compounds were determined from triplicate spikes of the individual compounds into separate buffer solutions.

Data from Reference 12.

TABLE 17

RECOVERY DATA FROM THREE LABORATORIES FOR SOLID-PHASE EXTRACTION (METHOD 3535)  
OF BASE/NEUTRAL/ACID EXTRACTABLES FROM SPIKED TCLP LEACHATES FROM SOIL SAMPLES

| Analyte               | Spike Level<br>( $\mu\text{g/L}$ )* | Lab 1 |     |    | Lab 2 |      |    | Lab 3 |     |    |
|-----------------------|-------------------------------------|-------|-----|----|-------|------|----|-------|-----|----|
|                       |                                     | %R    | RSD | n  | %R    | RSD  | n  | %R    | RSD | n  |
| o-Cresol              | 200,000                             | 86    | 8   | 7  | 35.3  | 0.7  | 3  | 7.6   | 6   | 3  |
| m-Cresol**            | --                                  | 77    | 8   | 7  | --    | --   | -- | --    | --  | -- |
| p-Cresol**            | --                                  | --    | --  | -- | --    | --   | -- | 7.7   | 11  | 3  |
| 2,4,6-Trichlorophenol | 2,000                               | 106   | 6   | 7  | 96.3  | 3.9  | 3  | 44.8  | 5   | 3  |
| 2,4,5-Trichlorophenol | 400,000                             | 93    | 3   | 7  | 80.5  | 4.5  | 3  | 63.3  | 11  | 3  |
| Pentachlorophenol     | 100,000                             | 79    | 2   | 7  | 33.8  | 12.2 | 3  | 29.2  | 13  | 3  |
| 1,4-Dichlorobenzene   | 7,500                               | 51    | 5   | 7  | 81.3  | 5.3  | 3  | 19.2  | 7   | 3  |
| Hexachloroethane      | 3,000                               | 50    | 5   | 7  | 66.2  | 2.1  | 3  | 12.6  | 11  | 3  |
| Nitrobenzene          | 2,000                               | 80    | 8   | 7  | 76.3  | 5.3  | 3  | 63.9  | 12  | 3  |
| Hexachlorobutadiene   | 500                                 | 53    | 8   | 7  | 63.3  | 4.8  | 3  | 9.6   | 9   | 3  |
| 2,4-Dinitrotoluene    | 130                                 | 89    | 8   | 7  | 35.7  | 2.6  | 3  | 58.2  | 17  | 3  |
| Hexachlorobenzene     | 130                                 | 84    | 21  | 7  | 92.3  | 1.6  | 3  | 71.7  | 9   | 3  |

(continued)

TABLE 17  
(continued)

| Analyte               | Spike Level (µg/L)* | Lab 1 |     |    | Lab 2 |      |    | Lab 3 |     |    |
|-----------------------|---------------------|-------|-----|----|-------|------|----|-------|-----|----|
|                       |                     | %R    | RSD | n  | %R    | RSD  | n  | %R    | RSD | n  |
| o-Cresol              | 200,00              | 97    | 13  | 7  | 37.8  | 4.5  | 3  | 6.1   | 24  | 3  |
| m-Cresol**            | --                  | 83    | 4   | 7  | --    | --   | -- | 6.0   | 25  | 3  |
| p-Cresol**            | --                  | --    | --  | -- | --    | --   | -- | --    | --  | -- |
| 2,4,6-Trichlorophenol | 2,000               | 104   | 4   | 7  | 91.7  | 8.0  | 3  | 37.7  | 25  | 3  |
| 2,4,5-Trichlorophenol | 400,000             | 94    | 4   | 7  | 85.2  | 0.4  | 3  | 64.4  | 10  | 3  |
| Pentachlorophenol     | 100,000             | 109   | 11  | 7  | 41.9  | 28.2 | 3  | 36.6  | 32  | 3  |
| 1,4-Dichlorobenzene   | 7,500               | 50    | 5   | 7  | 79.7  | 1.0  | 3  | 26.5  | 68  | 3  |
| Hexachloroethane      | 3,000               | 51    | 3   | 7  | 64.9  | 2.0  | 3  | 20.3  | 90  | 3  |
| Nitrobenzene          | 2,000               | 80    | 4   | 7  | 79.0  | 2.3  | 3  | 59.4  | 6   | 3  |
| Hexachlorobutadiene   | 500                 | 57    | 5   | 7  | 60    | 3.3  | 3  | 16.6  | 107 | 3  |
| 2,4-Dinitrotoluene    | 130                 | 86    | 6   | 7  | 38.5  | 5.2  | 3  | 62.2  | 6   | 3  |
| Hexachlorobenzene     | 130                 | 86    | 7   | 7  | 91.3  | 0.9  | 3  | 75.5  | 5   | 3  |

\* 250-mL aliquots of leachate were spiked. Lab 1 spiked at one-half these levels.

\*\* m-Cresol and p-Cresol coelute. Lab 1 and Lab 3 reported o-Cresol and the sum of — and p-Cresol. Lab 2 reported the sum of all three isomers of Cresol.

Data from Reference 12.

TABLE 18

SINGLE-LABORATORY PAH ANALYSIS DATA FROM A REAL SOIL CONTAMINATED WITH  
CREOSOTE, USING METHOD 3546  
(MICROWAVE EXTRACTION)

| Compound                | Concentration (µg/kg) | RSD (%) | REAC values (µg/kg) |
|-------------------------|-----------------------|---------|---------------------|
| Naphthalene             | 2,170                 | 12.4    | 710,000             |
| 2-Methylnaphthalene     | 28,710                | 3.1     | N/R                 |
| 1-Methylnaphthalene     | 33,180                | 2.4     | N/R                 |
| Biphenyl                | 13,440                | 6.0     | N/R                 |
| 2,6-Dimethylnaphthalene | 52,990                | 3.8     | N/R                 |
| Acenaphthylene          | 16,320                | 3.1     | 21,000              |
| Acenaphthene            | 801,210               | 6.0     | 1,700,000           |
| Fluorene                | 789,980               | 3.4     | 990,000             |
| Phenanthrene            | 1,627,480             | 0.7     | 3,300,000           |
| Anthracene              | 346,010               | 4.0     | 360,000             |
| Benzo(a)anthracene      | 300,380               | 2.7     | 310,000             |
| Fluoranthene            | 1,331,690             | 1.6     | 1,600,000           |
| Pyrene                  | 1,037,710             | 3.0     | 1,100,000           |
| Chrysene                | 293,200               | 3.4     | 320,000             |
| Benzo(b)fluoranthene    | 152,000               | 3.8     | 140,000             |
| Benzo(k)fluoranthene    | 127,740               | 3.6     | 130,000             |
| Benzo(e)pyrene          | 87,610                | 3.9     | N/R                 |
| Benzo(a)pyrene          | 128,330               | 3.9     | 110,000             |
| Perylene                | 35,260                | 4.3     | N/R                 |
| Indeno(123-cd)pyrene    | 63,900                | 5.0     | 25,000              |
| Dibenz(a,h)anthracene   | 17,290                | 6.9     | N/R                 |
| Benzo(ghi)perylene      | 42,720                | 6.9     | 20,000              |

\*n = 4

Soil samples obtained from US EPA Emergency Response Center archive bank through their contract laboratory REAC (Edison, NJ). The standard Soxhlet extraction procedures were performed by REAC three years earlier; this long storage period is believed to account for the low naphthalene recovery data in the present study

REAC data labeled N/R = not reported

TABLE 19

## SINGLE-LABORATORY PAH RECOVERY DATA FROM HS-5 MARINE SEDIMENT MATERIALS, USING METHOD 3546 (MICROWAVE EXTRACTION)

| Compound              | Certified Value (µg/kg) | Confidence Interval (µg/kg) | Recovery (%) |
|-----------------------|-------------------------|-----------------------------|--------------|
| Naphthalene           | 250                     | 180 - 320                   | 76           |
| Acenaphthylene        | 150                     | *                           | 107          |
| Acenaphthene          | 230                     | 130 - 330                   | 61           |
| Fluorene              | 400                     | 300 - 500                   | 63           |
| Phenanthrene          | 5,200                   | 4,200 - 6,200               | 72           |
| Anthracene            | 380                     | 230 - 530                   | 84           |
| Fluoranthene          | 8,400                   | 5,800 - 10,000              | 81           |
| Pyrene                | 5,800                   | 4,000 - 7,600               | 69           |
| Benzo(a)anthracene    | 2,900                   | 1,700 - 4,100               | 53           |
| Chrysene              | 2,800                   | 1,900 - 3,700               | 76           |
| Benzo(b)fluoranthene  | 2,000                   | 1,000 - 3,000               | 84           |
| Benzo(k)fluoranthene  | 1,000                   | 600 - 1,400                 | 137          |
| Benzo(a)pyrene        | 1,700                   | 900 - 2,500                 | 52           |
| Indeno(123-cd) pyrene | 1,300                   | 600 - 2,000                 | 63           |
| Dibenz(a,h)anthracene | 200                     | 100 - 300                   | 125          |
| Benzo(ghi)perylene    | 1,300                   | 1000 - 1600                 | 64           |

n = 3

\* values not certified

The uncertainties represent 90% confidence intervals

TABLE 20

## SINGLE-LABORATORY PAH RECOVERY DATA FROM HS-4 MARINE SEDIMENT MATERIALS, USING METHOD 3546 (MICROWAVE EXTRACTION)

| Compound              | Certified Value (µg/kg) | Confidence Interval (µg/kg) | Recovery (%) |
|-----------------------|-------------------------|-----------------------------|--------------|
| Naphthalene           | 150                     | *                           | 54           |
| Acenaphthylene        | 150                     | *                           | 82           |
| Acenaphthene          | 150                     | *                           | 63           |
| Fluorene              | 150                     | *                           | 81           |
| Phenanthrene          | 680                     | 600 - 760                   | 81           |
| Anthracene            | 140                     | 70 - 210                    | 108          |
| Fluoranthene          | 1250                    | 1,150 - 1,350               | 84           |
| Pyrene                | 940                     | 820 - 1,060                 | 85           |
| Benzo(a)anthracene    | 530                     | 470 - 580                   | 78           |
| Chrysene              | 650                     | 570 - 730                   | 84           |
| Benzo(b)fluoranthene  | 700                     | 550 - 850                   | 84           |
| Benzo(k)fluoranthene  | 360                     | 310 - 410                   | 156          |
| Benzo(a)pyrene        | 650                     | 570 - 730                   | 73           |
| Indeno(123-cd) pyrene | 510                     | 360 - 660                   | 88           |
| Dibenz(a,h)anthracene | 120                     | 70 - 170                    | 117          |
| Benzo(ghi)perylene    | 580                     | 360 - 800                   | 91           |

n = 3

\* values not certified

The uncertainties represent 90% confidence intervals

TABLE 21

## SINGLE-LABORATORY PAH RECOVERY DATA FROM HS-3 MARINE SEDIMENT MATERIALS, USING METHOD 3546 (MICROWAVE EXTRACTION)

| Compound              | Certified Value<br>( $\mu\text{g}/\text{kg}$ ) | Confidence Interval<br>( $\mu\text{g}/\text{kg}$ ) | Recovery (%) |
|-----------------------|--|--|--------------|
| Naphthalene           | 9,000  | 8300 - 9,700                                       | 61           |
| Acenaphthylene        | 300  | 200 - 400  | 199          |
| Acenaphthene          | 4,500  | 3,000 - 6,000                                      | 80           |
| Fluorene              | 13,300   | 10,200 -16,400                                     | 58           |
| Phenanthrene          | 85,000   | 65000 -105,000                                     | 87           |
| Anthracene            | 13,400   | 12,900 -13,900                                     | 48           |
| Fluoranthene          | 60,000   | 51,000-69,000                                      | 91           |
| Pyrene                | 39,000   | 30,000-48,000                                      | 86           |
| Benzo(a)anthracene    | 14,600   | 12,600-16,600                                      | 78           |
| Chrysene              | 14,100   | 12,100-16,100                                      | 91           |
| Benzo(b)fluoranthene  | 7,700  | 6,500-8,900  | 101          |
| Benzo(k)fluoranthene  | 2,800  | 800-4,800  | 275          |
| Benzo(a)pyrene        | 7,400  | 3,000-7,000  | 74           |
| Indeno(123-cd)pyrene  | 5,400  | 4,100-6,700  | 100          |
| Dibenz(a,h)anthracene | 1,300  | 800-1,800  | 118          |
| Benzo(ghi)perylene    | 5,000  | 3,000-7,000  | 99           |

n = 3

\* values not certified

The uncertainties represent 90% confidence intervals



TABLE 22


SINGLE-LABORATORY PAH RECOVERY DATA FROM SRM 1941 MARINE SEDIMENT,  
USING METHOD 3546 (MICROWAVE EXTRACTION)

| Compound              | Certified Value<br>(µg/kg) | Recovery<br>(%) |
|-----------------------|----------------------------|-----------------|
| Naphthalene           | 1010                       | 97.4            |
| Fluorene              | 100                        | 100.0           |
| Phenanthrene          | 490                        | 102.0           |
| Fluoranthene          | 980                        | 116.7           |
| Pyrene                | 810                        | 97.3            |
| Benz(a)anthracene     | 430                        | 89.8            |
| Chrysene              | 380                        | 130.3           |
| Benzo(b)fluoranthene  | 740                        | 95.8            |
| Benzo(k)fluoranthene  | 360                        | 130.2           |
| Benz(e)pyrene         | 550                        | 81.0            |
| Benzo(a)pyrene        | 630                        | 76.0            |
| Perylene              | 450                        | 72.4            |
| Indeno(123-cd)pyrene  | 500                        | 126.0           |
| Dibenz(a,h)anthracene | 110                        | 78.7            |
| Benz(ghi)perylene     | 530                        | 85.2            |

n = 3

All RSDs < 10%

FIGURE 1  
TAILING FACTOR CALCULATION




$$\text{TAILING FACTOR} = \frac{BC}{AB}$$

Example calculation: Peak Height = DE = 100 mm  
10% Peak Height = BD = 10 mm  
Peak Width at 10% Peak Height = AC = 23 mm  
AB = 11 mm  
BC = 12 mm

$$\text{Therefore: Tailing Factor} = \frac{12}{11} = 1.1$$

FIGURE 2  
GAS CHROMATOGRAM OF BASE/NEUTRAL AND ACID CALIBRATION STANDARD



## METHOD 6020A

### INDUCTIVELY COUPLED PLASMA-MASS SPECTROMETRY

SW-846 is not intended to be an analytical training manual. Therefore, method procedures are written based on the assumption that they will be performed by analysts who are formally trained in at least the basic principles of chemical analysis and in the use of the subject technology.

In addition, SW-846 methods, with the exception of required method use for the analysis of method-defined parameters, are intended to be guidance methods which contain general information on how to perform an analytical procedure or technique which a laboratory can use as a basic starting point for generating its own detailed Standard Operating Procedure (SOP), either for its own general use or for a specific project application. The performance data included in this method are for guidance purposes only, and are not intended to be and must not be used as absolute QC acceptance criteria for purposes of laboratory accreditation.

#### 1.0 SCOPE AND APPLICATION

1.1 Inductively coupled plasma-mass spectrometry (ICP-MS) is applicable to the determination of sub- $\mu\text{g/L}$  concentrations of a large number of elements in water samples and in waste extracts or digests (Refs. 1 and 2). When dissolved constituents are required, samples must be filtered and acid-preserved prior to analysis. No digestion is required prior to analysis for dissolved elements in water samples. Acid digestion prior to filtration and analysis is required for groundwater, aqueous samples, industrial wastes, soils, sludges, sediments, and other solid wastes for which total (acid-soluble) elements are required.

1.2 ICP-MS has been applied to the determination of over 60 elements in various matrices. Analytes for which EPA has demonstrated the acceptability of this method in a multi-laboratory study on solid and aqueous wastes are listed below.

| Element   |      | CASRN <sup>a</sup> |
|-----------|------|--------------------|
| Aluminum  | (Al) | 7429-90-5          |
| Antimony  | (Sb) | 7440-36-0          |
| Arsenic   | (As) | 7440-38-2          |
| Barium    | (Ba) | 7440-39-3          |
| Beryllium | (Be) | 7440-41-7          |
| Cadmium   | (Cd) | 7440-43-9          |
| Calcium   | (Ca) | 7440-70-2          |
| Chromium  | (Cr) | 7440-47-3          |
| Cobalt    | (Co) | 7440-48-4          |
| Copper    | (Cu) | 7440-50-8          |
| Iron      | (Fe) | 7439-89-6          |
| Lead      | (Pb) | 7439-92-1          |
| Magnesium | (Mg) | 7439-95-4          |
| Manganese | (Mn) | 7439-96-5          |
| Mercury   | (Hg) | 7439-97-6          |
| Nickel    | (Ni) | 7440-02-0          |

| Element   |      | CASRN <sup>a</sup> |
|-----------|------|--------------------|
| Potassium | (K)  | 7440-09-7          |
| Selenium  | (Se) | 7782-49-2          |
| Silver    | (Ag) | 7440-22-4          |
| Sodium    | (Na) | 7440-23-5          |
| Thallium  | (Tl) | 7440-28-0          |
| Vanadium  | (V)  | 7440-62-2          |
| Zinc      | (Zn) | 7440-66-6          |

<sup>a</sup>Chemical Abstract Service Registry Number

Acceptability of this method for an element was based upon the multi-laboratory performance compared with that of either furnace atomic absorption spectrophotometry or inductively coupled plasma-atomic emission spectrometry. It should be noted that one multi-laboratory study was conducted in 1988 and advances in ICP-MS instrumentation and software have been made since that time and additional studies have been added with validation and improvements in performance of the method. Performance, in general, exceeds the multi-laboratory performance data for the listed elements. It is expected that current performance will exceed the multi-laboratory performance data for the listed elements (and others) that are provided in Sec. 13.0. The lower limit of quantitation and linear ranges will vary with the matrices, instrumentation, and operating conditions. In relatively simple matrices, quantitation limits will generally be below 0.1 µg/L. Less sensitive elements (like Se and As) and desensitized major elements may be 1.0 µg/L or higher.

1.3 If this method is used to determine any analyte not listed in Sec. 1.2, it is the responsibility of the analyst to demonstrate the accuracy and precision of the method in the waste to be analyzed. The analyst is always required to monitor potential sources of interferences and take appropriate action to ensure data of known quality (see Sec. 9.0). Other elements and matrices may be analyzed by this method if performance is demonstrated for the analyte of interest, in the matrices of interest, at the concentration levels of interest in the same manner as the listed elements and matrices (see Sec. 9.0).

1.4 An appropriate internal standard is required for each analyte determined by ICP-MS. Recommended internal standards are <sup>6</sup>Li, <sup>45</sup>Sc, <sup>89</sup>Y, <sup>103</sup>Rh, <sup>115</sup>In, <sup>159</sup>Tb, <sup>165</sup>Ho, <sup>74</sup>Ge, and <sup>209</sup>Bi. The lithium internal standard should have an enriched abundance of <sup>6</sup>Li, so that interference from lithium native to the sample is minimized. Other elements may need to be used as internal standards when samples contain significant native amounts of the recommended internal standards.

1.5 Prior to employing this method, analysts are advised to consult the each preparative method that may be employed in the overall analysis (e.g., a 3000 series method) for additional information on quality control procedures, development of QC acceptance criteria, calculations, and general guidance. Analysts also should consult the disclaimer statement at the front of the manual and the information in Chapter Two for guidance on the intended flexibility in the choice of methods, apparatus, materials, reagents, and supplies, and on the responsibilities of the analyst for demonstrating that the techniques employed are appropriate for the analytes of interest, in the matrix of interest, and at the levels of concern.

In addition, analysts and data users are advised that, except where explicitly specified in a regulation, the use of SW-846 methods is *not* mandatory in response to Federal testing requirements. The information contained in this method is provided by EPA as guidance to be used by the analyst and the regulated community in making judgments necessary to generate results that meet the data quality objectives for the intended application.

1.6 Use of this method is restricted to use by, or under supervision of, properly experienced and trained personnel, including spectroscopists who are knowledgeable in the recognition and in the correction of spectral, chemical, and physical interferences in ICP-MS. Each analyst must demonstrate the ability to generate acceptable results with this method.

## 2.0 SUMMARY OF METHOD

2.1 Prior to analysis, samples should be solubilized or digested using the appropriate sample preparation methods (see Chapter Three). When analyzing groundwater or other aqueous samples for dissolved constituents, acid digestion is not necessary if the samples are filtered and acid preserved prior to analysis (refer to Sec. 1.1).

2.2 This method describes the multi-elemental determination of analytes by ICP-MS in environmental samples. The method measures ions produced by a radio-frequency inductively coupled plasma. Analyte species originating in a liquid are nebulized and the resulting aerosol is transported by argon gas into the plasma torch. The ions produced by high temperatures are entrained in the plasma gas and extracted through a differentially pumped vacuum interface and separated on the basis of their mass-to-charge ratio by a mass spectrometer. The ions transmitted through the mass spectrometer are quantified by a channel electron multiplier or Faraday detector and the ion information is processed by the instrument's data handling system. Interferences must be assessed and valid corrections applied or the data qualified to indicate problems. Interference correction must include compensation for background ions contributed by the plasma gas, reagents, and constituents of the sample matrix.

## 3.0 DEFINITIONS

Refer to Chapter One, Chapter Three, and the manufacturer's instructions for definitions that may be applicable to this procedure.

## 4.0 INTERFERENCES

4.1 Solvents, reagents, glassware, and other sample processing hardware may yield artifacts and/or interferences to sample analysis. All these materials must be demonstrated to be free from interferences under the conditions of the analysis by analyzing method blanks. Specific selection of reagents and purification of solvents by distillation in all-glass systems may be necessary. Refer to each method to be used for specific guidance on quality control procedures and to Chapter Three for general guidance on the cleaning of glassware. Also refer to the preparative methods to be used for discussions on interferences.

4.2 Isobaric elemental interferences in ICP-MS are caused by isotopes of different elements forming atomic ions with the same nominal mass-to-charge ratio ( $m/z$ ). A data system must be used to correct for these interferences. This involves determining the signal for another isotope of the interfering element and subtracting the appropriate signal from the analyte isotope signal. Since commercial ICP-MS instruments nominally provide unit resolution at 10% of the peak height, very high ion currents at adjacent masses can also contribute to ion signals at the mass of interest. Although this type of interference is uncommon, it is not easily corrected, and samples exhibiting a significant problem of this type could require resolution improvement, matrix separation, or analysis using another verified and documented isotope, or use of another method.

4.3 Isobaric molecular and doubly-charged ion interferences in ICP-MS are caused by ions consisting of more than one atom or charge, respectively. Most isobaric interferences that could affect ICP-MS determinations have been identified in the literature (Refs. 3 and 4). Examples include  $^{75}\text{ArCl}^+$  ion on the  $^{75}\text{As}$  signal and  $\text{MoO}^+$  ions on the cadmium isotopes. While the approach used to correct for molecular isobaric interferences is demonstrated below using the natural isotope abundances from the literature (Ref. 5), the most precise coefficients for an instrument can be determined from the ratio of the net isotope signals observed for a standard solution at a concentration providing suitable (<1%) counting statistics. Because the  $^{35}\text{Cl}$  natural abundance of 75.77% is 3.13 times the  $^{37}\text{Cl}$  abundance of 24.23%, the chloride correction for arsenic can be calculated (approximately) as follows (where the  $^{38}\text{Ar}^{37}\text{Cl}^+$  contribution at m/z 75 is a negligible 0.06% of the  $^{40}\text{Ar}^{35}\text{Cl}^+$  signal):

Corrected arsenic signal (using natural isotopes abundances for coefficient approximations) =  
(m/z 75 signal) - (3.13) (m/z 77 signal) + (2.73) (m/z 82 signal),

where the final term adjusts for any selenium contribution at 77 m/z.

NOTE: Arsenic values can be biased high by this type of equation when the net signal at m/z 82 is caused by ions other than  $^{82}\text{Se}^+$ , (e.g.,  $^{81}\text{BrH}^+$  from bromine wastes [Ref. 6]).

Similarly:

Corrected cadmium signal (using natural isotopes abundances for coefficient approximations) =  
(m/z 114 signal) - (0.027)(m/z 118 signal) - (1.63)(m/z 108 signal),

where last 2 terms adjust for any  $^{114}\text{Sn}^+$  or  $^{114}\text{MoO}^+$  contributions at m/z 114.

NOTE: Cadmium values will be biased low by this type of equation when  $^{92}\text{ZrO}^+$  ions contribute at m/z 108, but use of m/z 111 for Cd is even subject to direct ( $^{94}\text{ZrOH}^+$ ) and indirect ( $^{90}\text{ZrO}^+$ ) additive interferences when Zr is present.

NOTE: As for the arsenic equation above, the coefficients could be improved. The most appropriate coefficients for a particular instrument can be determined from the ratio of the net isotope signals observed for a standard solution at a concentration providing suitable (<1%) counting precision.

The accuracy of these types of equations is based upon the constancy of the observed isotopic ratios for the interfering species. Corrections that presume a constant fraction of a molecular ion relative to the "parent" ion have not been found (Ref. 7) to be reliable, e.g., oxide levels can vary with operating conditions. If a correction for an oxide ion is based upon the ratio of parent-to-oxide ion intensities, the correction must be adjusted for the degree of oxide formation by the use of an appropriate oxide internal standard previously demonstrated to form a similar level of oxide as the interferent. For example, this type of correction has been reported (Ref. 7) for oxide-ion corrections using  $\text{ThO}^+/\text{Th}^+$  for the determination of rare earth elements. The use of aerosol desolvation and/or mixed gas plasmas have been shown to greatly reduce molecular interferences (Ref. 8). These techniques can be used provided that the lower limits of quantitation, accuracy, and precision requirements for analysis of the samples can be met.

4.4 Additionally, solid phase chelation may be used to eliminate isobaric interferences from both element and molecular sources. An on-line method has been demonstrated for environmental waters such as sea water, drinking water and acid decomposed samples. Acid decomposed samples refer to samples decomposed by methods similar to Methods 3052, 3051, 3050 or 3015. Samples with percent levels of iron and aluminum should be avoided. The

method also provides a method for preconcentration to enhance quantitation limits simultaneously with elimination of isobaric interferences. The method relies on chelating resins such as imminodiacetate or other appropriate resins and selectively concentrates the elements of interest while eliminating interfering elements from the sample matrix. By eliminating the elements that are direct isobaric interferences or those that form isobaric interfering molecular masses, the mass region is simplified and these interferences can not occur. The method has been proven effective for the certification of standard reference materials and validated using SRMs (Refs. 13 through 15). The method has the potential to be used on-line or off-line as an effective sample preparation method specifically designed to address interference problems.

4.5 Physical interferences are associated with the sample nebulization and transport processes as well as with ion-transmission efficiencies. Nebulization and transport processes can be affected if a matrix component causes a change in surface tension or viscosity. Changes in matrix composition can cause significant signal suppression or enhancement (Ref. 9). Dissolved solids can deposit on the nebulizer tip of a pneumatic nebulizer and on the interface skimmers (reducing the orifice size and the instrument performance). Total solid levels below 0.2% (2,000 mg/L) are recommended (Ref. 10) to minimize solid deposition. An internal standard can be used to correct for physical interferences, if it is carefully matched to the analyte so that the two elements are similarly affected by matrix changes (Ref. 11). When intolerable physical interferences are present in a sample, a significant suppression of the internal standard signals (to less than 30% of the signals in the calibrations standard) will be observed. Dilution of the sample fivefold (1+4) will usually eliminate the problem (see Sec. 9.5).

4.6 Memory interferences or carry-over can occur when there are large concentration differences between samples or standards which are analyzed sequentially. Sample deposition on the sampler and skimmer cones, spray chamber design, and the type of nebulizer affect the extent of observed memory interferences. The rinse period between samples must be long enough to eliminate significant memory interference.

## 5.0 SAFETY

5.1 This method does not address all safety issues associated with its use. The laboratory is responsible for maintaining a safe work environment and a current awareness file of OSHA regulations regarding the safe handling of the chemicals specified in this method. A reference file of material safety data sheets (MSDSs) should be available to all personnel involved in these analyses.

5.2 Concentrated nitric and hydrochloric acids are moderately toxic and extremely irritating to skin and mucus membranes. Use these reagents in a hood and if eye or skin contact occurs, flush with large volumes of water. Always wear safety glasses or a shield for eye protection when working with these reagents. Hydrofluoric acid is a very toxic acid and penetrates the skin and tissues deeply if not treated immediately. Injury occurs in two stages; first, by hydration that induces tissue necrosis and then by penetration of fluoride ions deep into the tissue and by reaction with calcium. Boric acid and other complexing reagents and appropriate treatment agents should be administered immediately. Consult appropriate safety literature and have the appropriate treatment materials readily available prior to working with this acid. See Method 3052 for specific suggestions for handling hydrofluoric acid from a safety and an instrument standpoint.

5.3 Many metal salts are extremely toxic if inhaled or swallowed. Extreme care must be taken to ensure that samples and standards are handled properly and that all exhaust gases are properly vented. Wash hands thoroughly after handling.



5.4 The acidification of samples containing reactive materials may result in the release of toxic gases, such as cyanides or sulfides. For this reason, the acidification and digestion of samples should be performed in an approved fume hood.

## 6.0 EQUIPMENT AND SUPPLIES

6.1 Inductively coupled plasma-mass spectrometer -- A system capable of providing resolution, better than or equal to 1.0 amu at 10% peak height is required. The system must have a mass range from at least 6 to 240 amu and a data system that allows corrections for isobaric interferences and the application of the internal standard technique. Use of a mass-flow controller for the nebulizer argon and a peristaltic pump for the sample solution is recommended.

6.2 Argon gas supply -- High-purity grade (99.99%).

## 7.0 REAGENTS AND STANDARDS

7.1 Reagent- or trace metals-grade chemicals must be used in all tests. Unless otherwise indicated, it is intended that all reagents conform to the specifications of the Committee on Analytical Reagents of the American Chemical Society, where such specifications are available. Other grades may be used, provided it is first ascertained that the reagent is of sufficiently high purity to permit its use without lessening the accuracy of the determination.

7.2 Acids used in the preparation of standards and for sample processing must be of high purity. Redistilled acids are recommended because of the high sensitivity of ICP-MS. Nitric acid at less than 2% (v/v) is required for ICP-MS to minimize damage to the interface and to minimize isobaric molecular-ion interferences with the analytes. Many more molecular-ion interferences are observed when hydrochloric and sulfuric acids are used (Refs. 3 and 4). Concentrations of antimony and silver between 50-500 µg/L require 1% (v/v) HCl for stability; for concentrations above 500 µg/L Ag, additional HCl will be needed. Consequently, accuracy of analytes requiring significant chloride molecular ion corrections (such as As and V) will degrade.

7.3 Reagent water -- All references to water in the method refer to reagent water, unless otherwise specified. Reagent water must be free of interferences.

7.4 Standard stock solutions for each analyte may be purchased or prepared from ultra-high purity grade chemicals or metals (99.99 or greater purity). See Method 6010 for instructions on preparing standard solutions from solids.

7.4.1 Bismuth internal standard stock solution (1 mL = 100 µg of Bi) -- Dissolve 0.1115 g of Bi<sub>2</sub>O<sub>3</sub> in a minimum amount of dilute HNO<sub>3</sub>. Add 10 mL of conc. HNO<sub>3</sub> and dilute to 1,000 mL with reagent water.

7.4.2 Germanium internal standard stock solution (1 mL = 100 µg of Ge) -- Dissolve 0.2954 g of GeCl<sub>4</sub> in a minimum amount of dilute HNO<sub>3</sub>. Add 10 mL of conc. HNO<sub>3</sub> and dilute to 1,000 mL with reagent water.

7.4.3 Holmium internal standard stock solution (1 mL = 100 µg of Ho) -- Dissolve 0.1757 g of Ho<sub>2</sub>(CO<sub>3</sub>)<sub>2</sub>•5H<sub>2</sub>O in 10 mL of reagent water and 10 mL of HNO<sub>3</sub>. After dissolution is complete, warm the solution to degas. Add 10 mL conc. of HNO<sub>3</sub> and dilute to 1,000 mL with reagent water.

7.4.4 Indium internal standard stock solution (1 mL = 100 µg of In) -- Dissolve 0.1000 g of indium metal in 10 mL of conc. HNO<sub>3</sub>. Dilute to 1,000 mL with reagent water.

7.4.5 Lithium internal standard stock solution (1 mL = 100 µg of <sup>6</sup>Li) -- Dissolve 0.6312 g of 95-atom-% <sup>6</sup>Li, Li<sub>2</sub>CO<sub>3</sub> in 10 mL of reagent water and 10 mL of HNO<sub>3</sub>. After dissolution is complete, warm the solution to degas. Add 10 mL conc. of HNO<sub>3</sub> and dilute to 1,000 mL with reagent water.

7.4.6 Rhodium internal standard stock solution (1 mL = 100 µg of Rh) -- Dissolve 0.3593 g of ammonium hexachlororhodate (III) (NH<sub>4</sub>)<sub>3</sub>RhCl<sub>6</sub> in 10 mL reagent water. Add 100 mL of conc. HCl and dilute to 1,000 mL with reagent water.

7.4.7 Scandium internal standard stock solution (1 mL = 100 µg of Sc) -- Dissolve 0.15343 g of Sc<sub>2</sub>O<sub>3</sub> in 10 mL (1+1) of hot HNO<sub>3</sub>. Add 5 mL of conc. HNO<sub>3</sub> and dilute to 1,000 mL with reagent water.

7.4.8 Terbium internal standard stock solution (1 mL = 100 µg of Tb) -- Dissolve 0.1828 g of Tb<sub>2</sub>(CO<sub>3</sub>)<sub>3</sub>•5H<sub>2</sub>O in 10 mL (1+1) of HNO<sub>3</sub>. After dissolution is complete, warm the solution to degas. Add 5 mL of conc. HNO<sub>3</sub> and dilute to 1,000 mL with reagent water.

7.4.9 Yttrium internal standard stock solution (1 mL = 100 µg of Y) -- Dissolve 0.2316 g of Y<sub>2</sub>(CO<sub>3</sub>)<sub>3</sub>•3H<sub>2</sub>O in 10 mL (1+1) of HNO<sub>3</sub>. Add 5 mL conc. of HNO<sub>3</sub> and dilute to 1,000 mL with reagent water.

7.4.10 Titanium interference stock solution (1 mL = 100 µg of Ti) -- Dissolve 0.4133 g of (NH<sub>4</sub>)<sub>2</sub>TiF<sub>6</sub> in reagent water. Add 2 drops of conc. HF and dilute to 1,000 mL with reagent water.

7.4.11 Molybdenum interference stock solution (1 mL = 100 µg of Mo) -- Dissolve 0.2043 g of (NH<sub>4</sub>)<sub>2</sub>MoO<sub>4</sub> in reagent water. Dilute to 1,000 mL with reagent water.

7.4.12 Gold preservative stock solution for mercury (1 mL = 100 µg) -- Recommend purchasing as high purity prepared solution of AuCl<sub>3</sub> in dilute hydrochloric acid matrix.

7.5 Mixed calibration standard solutions are prepared by diluting the stock-standard solutions to levels in the linear range for the instrument in a solvent consisting of 1% (v/v) HNO<sub>3</sub> in reagent water. The calibration standard solutions must contain a suitable concentration of an appropriate internal standard for each analyte. Internal standards may be added on-line at the time of analysis using a second channel of the peristaltic pump and an appropriate mixing manifold. Generally, an internal standard should be no more than 50 amu removed from the analyte. Recommended internal standards include <sup>6</sup>Li, <sup>45</sup>Sc, <sup>89</sup>Y, <sup>103</sup>Rh, <sup>115</sup>In, <sup>159</sup>Tb, <sup>169</sup>Ho, <sup>74</sup>Ge and <sup>209</sup>Bi. Prior to preparing the mixed standards, each stock solution must be analyzed separately to determine possible spectral interferences or the presence of impurities. Care must be taken when preparing the mixed standards to ensure that the elements are compatible and stable together. Transfer the mixed standard solutions to freshly acid-cleaned FEP fluorocarbon or previously unused polyethylene or polypropylene bottles for storage. For all intermediate and working standards, especially low level standards (i.e., <1 ppm), stability must be demonstrated prior to use. Fresh mixed standards must be prepared as needed with the realization that concentrations can change on aging. (Refer to Sec. 10.3.1 for guidance on determining the viability of standards.)

7.6 Blanks -- Three types of blanks are required for the analysis. The calibration blank is used in establishing the calibration curve. The method blank is used to monitor for

possible contamination resulting from either the reagents (acids) or the equipment used during sample processing including filtration. The rinse blank is used to flush the system between all samples and standards.

7.6.1 The calibration blank consists of the same concentration(s) of the same acid(s) used to prepare the final dilution of the calibrating solutions of the analytes [often 1% HNO<sub>3</sub> (v/v) in reagent water] along with the selected concentrations of internal standards such that there is an appropriate internal standard element for each of the analytes. Use of HCl for antimony and silver is cited in Sec. 7.2.

7.6.2 The method blank must contain all of the reagents in the same volumes as used in the processing of the samples. The method blank must be carried through the complete procedure and contain the same acid concentration in the final solution as the sample solution used for analysis (refer to Sec. 9.9).

7.6.3 The rinse blank consists of 1 to 2% of HNO<sub>3</sub> (v/v) in reagent water. Prepare a sufficient quantity to flush the system between standards and samples. If mercury is to be analyzed, the rinse blank should also contain 2 µg/mL (ppm) of AuCl<sub>3</sub> solution.

7.7 The interference check solution (ICS) is prepared to contain known concentrations of interfering elements that will demonstrate the magnitude of interferences and provide an adequate test of any corrections. Chloride in the ICS provides a means to evaluate software corrections for chloride-related interferences such as <sup>35</sup>Cl<sup>16</sup>O<sup>+</sup> on <sup>51</sup>V<sup>+</sup> and <sup>40</sup>Ar<sup>35</sup>Cl<sup>+</sup> on <sup>75</sup>As<sup>+</sup>. Iron is used to demonstrate adequate resolution of the spectrometer for the determination of manganese. Molybdenum serves to indicate oxide effects on cadmium isotopes. The other components are present to evaluate the ability of the measurement system to correct for various molecular-ion isobaric interferences. The ICS is used to verify that the interference levels are corrected by the data system within quality control limits.

NOTE: The final ICS solution concentrations in Table 1 are intended to evaluate corrections for known interferences on only the analytes in Sec. 1.2. If this method is used to determine an element not listed in Sec. 1.2, the analyst should modify the ICS solutions, or prepare an alternative ICS solution, to allow adequate verification of correction of interferences on the unlisted element (see Sec. 9.7).

7.7.1 These solutions must be prepared from ultra-pure reagents. They can be obtained commercially or prepared by the following procedure.

7.7.1.1 Mixed ICS solution I may be prepared by adding 13.903 g of Al(NO<sub>3</sub>)<sub>3</sub>•9H<sub>2</sub>O, 2.498 g of CaCO<sub>3</sub> (dried at 180 °C for 1 hr before weighing), 1.000 g of Fe, 1.658 g of MgO, 2.305 g of Na<sub>2</sub>CO<sub>3</sub>, and 1.767 g of K<sub>2</sub>CO<sub>3</sub> to 25 mL of reagent water. Slowly add 40 mL of (1+1) HNO<sub>3</sub>. After dissolution is complete, warm the solution to degas. Cool and dilute to 1,000 mL with reagent water.

7.7.1.2 Mixed ICS solution II may be prepared by slowly adding 7.444 g of 85 % H<sub>3</sub>PO<sub>4</sub>, 6.373 g of 96% H<sub>2</sub>SO<sub>4</sub>, 40.024 g of 37% HCl, and 10.664 g of citric acid C<sub>6</sub>O<sub>7</sub>H<sub>8</sub> to 100 mL of reagent water. Dilute to 1,000 mL with reagent water.

7.7.1.3 Mixed ICS solution III may be prepared by adding 1.00 mL each of 100-µg/mL arsenic, cadmium, selenium, chromium, cobalt, copper, manganese, nickel, silver, vanadium, and zinc stock solutions to about 50 mL of

reagent water. Add 2.0 mL of concentrated HNO<sub>3</sub>, and dilute to 100.0 mL with reagent water.

#### 7.7.1.4 Working ICS solutions

7.7.1.4.1 ICS-A may be prepared by adding 10.0 mL of mixed ICS solution I (Sec. 7.7.1.1), 2.0 mL each of 100-µg/mL titanium stock solution (Sec. 7.4.9) and molybdenum stock solution (Sec. 7.4.10), and 5.0 mL of mixed ICS solution II (Sec. 7.7.1.2). Dilute to 100 mL with reagent water. ICS solution A must be prepared fresh weekly.

7.7.1.4.2 ICS-AB may be prepared by adding 10.0 mL of mixed ICS solution I (Sec. 7.7.1.1), 2.0 mL each of 100-µg/mL titanium stock solution (Sec. 7.4.9) and molybdenum stock solution (Sec. 7.4.10), 5.0 mL of mixed ICS solution II (Sec. 7.7.1.2), and 2.0 mL of mixed ICS solution III (Sec. 7.7.1.3). Dilute to 100 mL with reagent water. Although the ICS solution AB must be prepared fresh weekly, the analyst should be aware that the solution may precipitate silver more quickly.

7.8 The initial calibration verification (ICV) standard is prepared by the analyst (or a purchased second source reference material) by combining compatible elements from a standard source different from that of the calibration standard, and at concentration near the midpoint of the calibration curve (see Sec. 10.4.3 for use). This standard may also be purchased.

7.9 The continuing calibration verification (CCV) standard should be prepared in the same acid matrix using the same standards used for calibration, at a concentration near the mid-point of the calibration curve (see Sec. 10.4.4 for use).

7.10 Mass spectrometer tuning solution. A solution containing elements representing all of the mass regions of interest (for example, 10 µg/L of Li, Co, In, and Tl) must be prepared to verify that the resolution and mass calibration of the instrument are within the required specifications (see Sec. 10.2). This solution is also used to verify that the instrument has reached thermal stability (see Sec. 11.4).

## 8.0 SAMPLE COLLECTION, PRESERVATION, AND STORAGE

8.1 See the introductory material in Chapter Three, "Inorganic Analytes."

8.2 Only polyethylene or fluorocarbon (TFE or PFA) containers are recommended for use in this method.

## 9.0 QUALITY CONTROL

9.1 Refer to Chapter One for additional guidance on quality assurance (QA) and quality control (QC) protocols. When inconsistencies exist between QC guidelines, method-specific QC criteria take precedence over both technique-specific criteria and those criteria given in Chapter One, and technique-specific QC criteria take precedence over the criteria in Chapter One. Any effort involving the collection of analytical data should include development of a structured and systematic planning document, such as a Quality Assurance Project Plan (QAPP) or a Sampling and Analysis Plan (SAP), which translates project objectives and specifications into directions for those that will implement the project and assess the results.

Each laboratory should maintain a formal quality assurance program. The laboratory should also maintain records to document the quality of the data generated. All data sheets and quality control data should be maintained for reference or inspection.

9.2 Refer to a 3000 series method (Method 3005, 3010, 3015, 3031, 3040, 3050, 3051, or 3052) for appropriate QC procedures to ensure the proper operation of the various sample preparation techniques.

9.3 Instrument detection limits (IDLs) are a useful tool to evaluate the instrument noise level and response changes over time for each analyte from a series of reagent blank analyses to obtain a calculated concentration. They are not to be confused with the lower limits of quantitation, nor should they be used in establishing this limit. It may be helpful to compare the calculated IDLs to the established lower limit of quantitation, however, it should be understood that the lower limit of quantitation needs to be verified according to the guidance in Sec. 10.2.3.

IDLs in  $\mu\text{g/L}$  can be estimated by calculating the average of the standard deviations of three runs on three non-consecutive days from the analysis of a reagent blank solution with seven consecutive measurements per day. Each measurement should be performed as though it were a separate analytical sample (i.e., each measurement must be followed by a rinse and/or any other procedure normally performed between the analysis of separate samples). IDLs should be determined at least every three months or at a project-specific designated frequency and kept with the instrument log book. Refer to Chapter One for additional guidance.

#### 9.4 Initial demonstration of proficiency

Each laboratory must demonstrate initial proficiency with each sample preparation (a 3000 series method) and determinative method combination it utilizes by generating data of acceptable accuracy and precision for target analytes in a clean matrix. If an autosampler is used to perform sample dilutions, before using the autosampler to dilute samples, the laboratory should satisfy itself that those dilutions are of equivalent or better accuracy than is achieved by an experienced analyst performing manual dilutions. The laboratory must also repeat the demonstration of proficiency whenever new staff members are trained or significant changes in instrumentation are made.

9.5 Dilute and reanalyze samples that exceed the linear dynamic range or use an alternate, less sensitive calibration for which quality control data are already established.

9.6 The intensities of all internal standards must be monitored for every analysis. If the intensity of any internal standard in a sample falls below 70% of the intensity of that internal standard in the initial calibration standard, a significant matrix effect must be suspected. As an example, if the initial calibration internal standard response is 100,000 cps, anything below 70,000 cps in the sample would be unacceptable. Under these conditions, the established lower limit of quantitation has degraded and the correction ability of the internal standardization technique becomes questionable. The following procedure is followed -- First, make sure the instrument has not drifted by observing the internal standard intensities in the nearest clean matrix (calibration blank, Sec. 7.6.1). If the low internal standard intensities are also seen in the nearest calibration blank, terminate the analysis, correct the problem, recalibrate, verify the new calibration, and reanalyze the affected samples. If drift has not occurred, matrix effects need to be removed by dilution of the affected sample. The sample must be diluted fivefold (1+4) and reanalyzed with the addition of appropriate amounts of internal standards. If the first dilution does not eliminate the problem, this procedure must be repeated until the internal-standard intensities rise to the minimum 70% limit. Reported results must be corrected for all dilutions.

9.7 To obtain analyte data of known quality, it is necessary to measure more than the analytes of interest in order to apply corrections or to determine whether interference corrections are necessary. For example, tungsten oxide moleculars can be very difficult to distinguish from mercury isotopes. If the concentrations of interference sources (such as C, Cl, Mo, Zr, W) are such that, at the correction factor, the analyte is less than the limit of quantification and the concentration of interferents are insignificant, then the data may go uncorrected. Note that monitoring the interference sources does not necessarily require monitoring the interferant itself, but that a molecular species may be monitored to indicate the presence of the interferent. When correction equations are used, all QC criteria must also be met. Extensive QC for interference corrections are required at all times. The monitored masses must include those elements whose hydrogen, oxygen, hydroxyl, chlorine, nitrogen, carbon and sulfur molecular ions could impact the analytes of interest. Unsuspected interferences may be detected by adding pure major matrix components to a sample to observe any impact on the analyte signals. When an interference source is present, the sample elements impacted must be flagged to indicate (a) the percentage interference correction applied to the data or (b) an uncorrected interference by virtue of the elemental equation used for quantitation. The isotope proportions for an element or molecular-ion cluster provide information useful for quality assurance.

NOTE: Only isobaric elemental, molecular, and doubly charged interference corrections which use the observed isotopic-response ratios or parent-to-oxide ratios (provided an oxide internal standard is used as described in Sec. 4.2) for each instrument system are acceptable corrections for use in Method 6020.

9.8 For each batch of samples processed, at least one method blank must be carried throughout the entire sample preparation and analytical process, as described in Chapter One. A method blank is prepared by using a volume or weight of reagent water at the volume or weight specified in the preparation method, and then carried through the appropriate steps of the analytical process. These steps may include, but are not limited to, prefiltering, digestion, dilution, filtering, and analysis. If the method blank does not contain target analytes at a level that interferes with the project-specific DQOs, then the method blank would be considered acceptable.

In the absence of project-specific DQOs, if the blank is less than 10% of the lower limit of quantitation check sample concentration, less than 10% of the regulatory limit, or less than 10% of the lowest sample concentration for each analyte in a given preparation batch, whichever is greater, then the method blank is considered acceptable. If the method blank cannot be considered acceptable, the method blank should be re-run once, and if still unacceptable, then all samples after the last acceptable method blank should be reprepared and reanalyzed along with the other appropriate batch QC samples. These blanks will be useful in determining if samples are being contaminated. If the method blank exceeds the criteria, but the samples are all either below the reporting level or below the applicable action level or other DQOs, then the sample data may be used despite the contamination of the method blank.

#### 9.9 Laboratory control sample (LCS)

For each batch of samples processed, at least one LCS must be carried throughout the entire sample preparation and analytical process. The laboratory control samples should be spiked with each analyte of interest at the project-specific action level or, when lacking project-specific action levels, at approximately mid-point of the linear dynamic range. Acceptance criteria should either be defined in the project-specific planning documents or set at a laboratory derived limit developed through the use of historical analyses. In the absence of project-specific or historical data generated criteria, this limit should be set at  $\pm 20\%$  of the spiked value. Acceptance limits derived from historical data should be no wider than  $\pm 20\%$ . If the laboratory control sample is not acceptable, then the laboratory control sample should be re-run once and,

if still unacceptable, all samples after the last acceptable laboratory control sample should be reprepared and reanalyzed.

Concurrent analyses of standard reference materials (SRMs) containing known amounts of analytes in the media of interest are recommended and may be used as an LCS. For solid SRMs, 80 - 120% accuracy may not be achievable and the manufacturer's established acceptance criterion should be used for soil SRMs.

#### 9.10 Matrix spike, unspiked duplicate, or matrix spike duplicate (MS/Dup or MS/MSD)

Documenting the effect of the matrix, for a given preparation batch consisting of similar sample characteristics, should include the analysis of at least one matrix spike and one duplicate unspiked sample or one matrix spike/matrix spike duplicate pair. The decision on whether to prepare and analyze duplicate samples or a matrix spike/matrix spike duplicate pair should be based on a knowledge of the samples in the sample batch or as noted in the project-specific planning documents. If samples are expected to contain target analytes, then laboratories may use one matrix spike and a duplicate analysis of an unspiked field sample. If samples are not expected to contain target analytes, laboratories should use a matrix spike and matrix spike duplicate pair.

For each batch of samples processed, at least one MS/Dup or MS/MSD sample set should be carried throughout the entire sample preparation and analytical process as described in Chapter One. MS/MSDs are intralaboratory split samples spiked with identical concentrations of each analyte of interest. The spiking occurs prior to sample preparation and analysis. An MS/Dup or MS/MSD is used to document the bias and precision of a method in a given sample matrix.

Refer to Chapter One for definitions of bias and precision, and for the proper data reduction protocols. MS/MSD samples should be spiked at the same level, and with the same spiking material, as the corresponding laboratory control sample that is at the project-specific action level or, when lacking project-specific action levels, at approximately mid-point of the linear dynamic range. Acceptance criteria should either be defined in the project-specific planning documents or set at a laboratory-derived limit developed through the use of historical analyses per matrix type analyzed. In the absence of project-specific or historical data generated criteria, these limits should be set at  $\pm 25\%$  of the spiked value for accuracy and 20 relative percent difference (RPD) for precision. Acceptance limits derived from historical data should be no wider than  $\pm 25\%$  for accuracy and 20% for precision. Refer to Chapter One for additional guidance. If the bias and precision indicators are outside the laboratory control limits, if the percent recovery is less than 75% or greater than 125%, or if the relative percent difference is greater than 20%, then the interference test discussed in Sec. 9.11 should be conducted.

9.10.1 The relative percent difference between spiked matrix duplicate or unspiked duplicate determinations is to be calculated as follows:

$$RPD = \frac{|D_1 - D_2|}{\left(\frac{D_1 + D_2}{2}\right)} \times 100$$

where:

RPD = relative percent difference.